



Operating Manual

Benchttop Compact Series Priorclave
Laboratory Autoclaves



Introduction

Priorclave autoclaves are a range of compact laboratory autoclaves intended primarily for preparation, the making safe of ordinary laboratory and pathological waste and sterilisation purposes. The autoclaves are manufactured to a high standard and feature a number of design features. The sophisticated microprocessor control system provides a very simple of setting even the most complex cycles. The machines have been designed from the outset for safe operation and maintenance.

Properly looked after and serviced your autoclave should give years of valuable and trouble

Priorclave Service

Model No: _____ *

SerialNumber: _____

NBBI No: _____

Date of Manufacture: _____

Software Version: _____

Please quote the above when asking for parts or service:

* Model Numbers are in the format PS/Mxx/znnn
Mxx MID-Non-Vacuum models / MA - Vacuum Models
z V for vertical, H for Horizontal indicating Compact model, either 40 litre horizontal or vertical. (Horizontal 60 litre units have this position)
nnn Nominal working volume in Litres

PRIORCLAVE LIMITED
 129-131 Nathan Way
 West Thames Road Business Park
 London
 SE28 OAB
 Telephone: +44 (0)203166620
 Fax: +44 (0)203550616
 E-mail: service@priorclave.co.uk
 Website: www.priorclave.co.uk

or your local service agent:

Contents

Introduction.....	5
Notices and Important Information.....	8
Symbols used on the product and in this manual and their meanings.....	8
General Safety Instructions.....	10
Specification Table.....	11
Important Notices and Warnings.....	12
Operating Summary.....	16
Cycle Abort and Thermal Lock Override.....	16
Layout Diagram.....	17
Operation.....	19
Opening the Door.....	19
Checking Water Level.....	19
Loading.....	20
Settings.....	21
Automatic Timed Free Steaming.....	22
Pulsed Free Steaming.....	22
Load Sensed Process Timing.....	23
Rapid Cooling.....	23
Media Warming.....	24
Delayed Start Time.....	24
Vacuum Options (MVA (van) models only).....	24
Multi Program Memory Options.....	25
Closing the pressure door.....	26
Starting a cycle.....	26
Vent button.....	27
During the process time.....	27
Power loss during a running cycle.....	28
Cooling.....	28
Thermal lock.....	28
Cycle complete.....	29
Media Warming.....	29
Aborting a Cycle.....	29
Emergency «E»Stop Button (if fitted).....	29
Operation with Options & Accessories.....	30
Setting Lock Key Switch Option.....	30
Printer.....	30
Serial Interface.....	31

USB Interface.....	31
Automatic Waterfill Option.....	31
Air Intake Filter.....	32
Vent Filter Option.....	32
Internal Validation System.....	32
Accelerated media cooling.....	32
Chart Recorder / Data Recorder / Data..Logger.....	32
Changing Date & Time.....	33
Warning Indicators and Fault Codes.....	34
Routine Operator Maintenance.....	36
Daily Maintenance.....	36
Weekly Maintenance.....	36
Monthly Maintenance.....	37
Bi-Annual Maintenance.....	38
Fault Finding & Rectification Guide.....	41
Steam Table.....	43
Notes.....	44
Other Items Fitted.....	45

Notices and Important Information

Symbols used on the product and in this manual and their meanings

WARNING: Mechanical Hazard

In this manual, warnings draw attention to the potential for Danger to person including risk of severe injury or death. Each Mechanical Hazard Warning is emphasized by this icon.



WARNING: Electrical Hazard

In this manual, warnings draw attention to the potential for Danger to person including risk of severe injury or death. Each Electrical Hazard Warning is emphasized by this icon.



WARNING: Bio-hazard

In this manual, warnings draw attention to the potential for Danger to person including risk of severe injury or death. Each Bio-hazard Warning is emphasized by this icon.



Caution: -Heavy

In this manual, cautions draw attention to the potential for injury to personnel where product or item has a weight of over 25kg. Reference should be made to the instructions for heavy lifting before attempting to move or lift. Each Caution Heavy Warning is emphasized by this icon.



Caution: -Please Note

In this manual, cautions draw attention to the potential for Damage to equipment. Each Caution Please Note Warning is emphasized by this icon.



Hazard-Isolate before Access

When one of these stickers has been placed on a removable panel the power must be turned off before the panel is removed. There may be a number of areas behind the label which constitute a hazard. All such panels are service access panels only and should not be removed unless there is a full understanding of the equipment.



Electrical Hazard-Isolate before Access

When one of these stickers has been placed on a removable panel the power must be turned off before the panel is removed. There may be a number of areas behind the label which constitute an electrical shock hazard. All such panels are service access panels only and should not be removed unless there is a full understanding of the equipment.



Electrical Earth Point

This protective label indicates a point at which an electrical earth cable should be connected. When removing and replacing panels after maintenance electrical earth cables should be reconnected at these points.



Mechanical Hazard

When one of these stickers has been placed on a removable panel the power must be turned off before the panel is removed. There may be a number of areas behind the label which constitute a mechanical hazard. All such panels are service access panels only and should not be removed unless there is a full understanding of the equipment.



Caution Hot Surface.

Parts of the surface to which one of these stickers are attached may become hot during the operation of the equipment. Take care not to touch these surfaces without heat protection.



Refer to Manual



General Safety Instructions

All cleaning and servicing requires the autoclave to be isolated from the power source and the gas supply. This equipment weighs in excess of 18 kg (40 lbs) and requires at least 2 persons to lift.

Should a fault occur with the autoclave, immediately isolate and disconnect the incoming gas supply. When the autoclave is being transported it should be sat firmly on a trolley or pallet and be strapped upright to a flat pallet during transport.

When the machine is not in use, where possible it should be disconnected from the electrical and gas water supply.

Ensure the equipment is installed, operated and maintained by trained personnel.

Always isolate the autoclave before cleaning or maintenance

Always ensure the machine is level when in use.

If fitted, in the event of an emergency press the stop button situated on the front of the autoclave.

The autoclave should be used as provided and should not be tampered with or altered as it contains inbuilt safety systems, which could be compromised by any interference.


During operation some autoclave surfaces may become uncomfortably hot. Take care if touching these surfaces without thermal protection.

Care should be taken when opening the autoclave after a sterilising cycle as it will be hot and steam may be released. Heatproof gloves and a face shield should always be worn when unloading.

When operating the autoclave contains steam at elevated temperature and pressure. Always follow the instructions when operating autoclaves.

The autoclave should only be used for its intended purpose. You must consult the manufacturer before using the autoclave for anything other than its intended purpose.

Specification Table

Model:	PS/MID/C40	PS/MVA/C40	PS/MID/H60	PS/MVA/H60	
Intended use.	Principle use: Preparation of Laboratory growth media and sterilisation of laboratory waste prior to disposal. Additional uses: Sterilisation of bottle liquids and glassware, sterilisation of unwrapped instruments porous loads. * dependent upon options fitted				
	This equipment has not been designed for or intended for use as a medical steriliser.				
Capacity (L)	40	40	60	60	
Dimensions	External (mm) (wxdxh)				
	520x610x640	695x610x640	520x815x640	695x815x640	
Weight (Approx.)	Unloaded	90	95	95	102
	With water charge	96	101	102	109
	Door	20	20	20	20
	Pressure Vessel	10	10	14	14
Electrical	Voltage(V)	230			
	Frequency (Hz)	50/60			
	Phases	Single Phase with earth			
	Max. Current Rating (A)	14			
	with Drying(A)	14			
	Max.Heater power kW	2.9			
Steam	Steam Supply Required	Not Applicable			
Water	Softened Water Supply Required	Hand Fill			
	With Auto water fill	15mm BSP			
Drainage	Main Drain Connection	15mm BSP			
	Secondary Drain Connections (where Applicable)	15mm BSP			
Max. Sound Level	Without Vacuum Pump	Approx. 30db @ 1M			
	With Vacuum Pump	Approx. 65db @ 1M			
Max. Heat Emission	Full Cycle to thermal lock temperature*	$(2.9 \times 1.45 \times p) \times 0.75 / t$ kW/Hour Where: h = heat time (hrs) p = process time (hrs) t = Total cycle* (hrs)			

#Must be fitted with option PC/MID/PS. Please refer your order paper work or the specifications sheet with this manual

Important Notices and Warnings



Before despatch from our works all Priorclaves are subjected to rigorous electric the appropriate standards. Should you or your contractors carry out further insu tests as part of your internal procedures please disconnect the switch mode power testing. Failure to do so will result in a test failure and may lead to corruption of microprocessor memory which cannot be covered by our warranty.

Safety

If you are unclear about any aspects of this manual, the use and operation of the your autoclave process please contact Priorclave or your authorised Priorclave de proceeding.



Always wear gloves, a facemask and adequate protective clothing when unloading an autoclave and ensure that the workload does not exceed safe limits.

Priorclave are pleased to arrange training for operators in the use of their autoclaves at an extra charge.

Thermal Lock

The safety Thermal Lock (door retention device) has been set in accordance with and procedure defined in paras. 3.3.3.2.3 and 3.3.3.3 of BS2646 Part 5:1993.



The relatively light load defined under this procedure may not be appropriate to the autoclave in your Priorclave. Therefore, to ensure compliance with Health & Safety Guidance Note PM73 'Safety at Autoclaves' you should only advise to have your autoclave with its normal working load formally validated, and the lift up accordingly by properly trained personnel.

Stainless Steel Pressure Vessels.

Vessels are manufactured from grade 316 stabilised stainless steel, designed built in accordance with PD5500 category 3 as required by BS2646 Part 1.



Grade 316 stainless steel is employed to reduce the corrosive effects of substances such as hydroxides and chlorine. However we recommend that the interior of the vessel is kept free of such potentially harmful substances and is regularly cleaned out with soft water. The use of chlorine based or other aggressive cleaners is not recommended. Exposure to such substances could damage the surface finish and the integrity of the pressure vessel and do not use such cleaners. It should also be taken not to routinely introduce such chemicals where they are used to clean items that form part of the load. In such cases the items should be cleaned before autoclaving.

Product Life

Due to fatigue occurring in normal use the life of all pressure vessels is finite due to corrosion, erosion or other damage. Using calculations from PD5500, and assuming the maximum working pressure of 2.4 bar this gives the autoclave a predicted fatigue life of 15,000 operating cycles. The lifespan of the autoclave will obviously depend on the frequency of use, but for example (based on a 365 day working year) if the autoclave is used two or four times per day this gives a working life of 15 to 10.2 years respectively. You should consider the usage of the autoclave to determine the actual lifespan of the autoclave.

Cleaning



External cleaning should only be carried out with a damp cloth or with proprietary non-abrasive cleaners.

Water Supply and Quality



This autoclave has been designed to operate most effectively with softened water. Unless this autoclave has been specifically adapted for purified water supplies then demineralised, distilled or RO water supplies MUST NOT be used as the controls fitted measure electrical conductivity to detect water levels.



Connection to a hard water supply can lead to a build up of scale and will damage the heating and other parts of the system and could invalidate the warranty.

RO and Ultra Pure water can also damage some elements of the steam generation system unless the autoclave has been specified and modified to operate with water of high purity. The usual method of filling is hot.

* Please refer to the specification sheet included with this manual for details.

Servicing and Maintenance of Priorclave Autoclaves



Priorclave Laboratory Autoclaves are complex pressure systems designed to special regulations and as such should only be serviced or maintained by properly trained personnel.

If your autoclave is run at an average frequency of more than 3 times per week then we strongly recommend that it should be serviced every six months, even during its initial 12 month manufacturer's warranty period to maintain it in peak operating condition. Autoclaves used less frequently can be serviced at 12 month intervals.

Service contracts for preventative routine maintenance can be arranged with Priorclave (service@priorclave.co.uk) or with your Priorclave authorised service agent.

Priorclave Ltd. cannot be held responsible for hazards or damage resulting from the use of the pressure system including its closure components by untrained or unauthorised personnel. If in doubt please contact Priorclave (service@priorclave.co.uk) or your nearest authorised service agent.



Faults caused by servicing by unauthorised service agents will not be covered by any warranty supplied with the autoclave.

CE Marking

The CE mark applied to this autoclave is applied in relation to the EMC (Electromagnetic Compatibility) directive and the Low Voltage directive of the European Community. This indicates that this Priorclave autoclave meets the following technical standards:

~~BSEN61000-3~~

Electromagnetic Compatibility. Generic Emission Standard. Residential, Commercial and Industrial.

~~BSEN61000-1~~

Electromagnetic Compatibility. Generic Immunity Standard. Residential, Commercial and Industrial.

~~BSEN61010~~

Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Processes.

~~BSEN61010-040~~

Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Processes. Part 41, Particular Requirements for Autoclaves using Steam for the treatment of Materials and for Laboratory Processes.

Pressure vessels

PD5500

Unfired fusion welded pressure vessels

BS2646 1996

Autoclaves for sterilisation in laboratories

Conformity assessment modules B1 + D of the European Pressure equipment directive applied to ensure compliance with essential safety requirements.

A "Declaration of Conformity" in accordance with the CE marking requirements is available as a separate file at:

Priorclave Ltd.
129 /131 Nathan Way
West Thamesmead Business Park
London SE28 0AB

Environmental Conditions

This equipment has been designed for safe operation within the following environmental conditions:

< Indoor Use.

< Altitude up to 2,000 - 6,500 feet

At this altitude the standard safety valve set at 15 bar will be the upper limit of the operating pressure. Adjustments can be made to compensate for the required pressure. See Appendix B Steam Table for special conditions affecting calibration operation at elevated altitudes.

< Operating temperatures between 15°C and 27°C. - The cooling performance of air cooled autoclaves, however will be significantly affected at higher temperature range.

< Maximum Relative Humidity of 85% at any temperature between 15°C and 27°C, Non Condensing

< Mains Supply Voltage Variations not exceeding 10% of that shown on the Serial

Electromagnetic Interference

This equipment has been designed to comply with the requirements for immunity to electromagnetic interference under normal conditions of use. Care should be taken when positioning the equipment however, to avoid interference from potential external sources such as MR scanners or other equipment.

Quick Opening Doors



Extracts from NONMANDATORY APPENDIX FF (GUIDE FOR THE DESIGN AND OPERATION OF QUICK-ACTUATING (QUICK-OPENING) CLOSURES) from ASME SECTION VIII DIVISION 1

FF-6 INSPECTION

It is recommended that the user inspect the completed installation including the gauges before it is permitted to operate. Records of this inspection should be retained.

It is recommended that the user establishes, and documents a periodic inspection program, and that, this program is followed and documented.

FF-7 TRAINING

Many accidents involving quick-closing closures have occurred because the operators have been unfamiliar with the equipment or its safety features. The greater safety of current designs has sometimes been produced by the use of sophisticated mechanical and electronic control devices. To ensure that these features produce the maximum safety, personnel should be properly trained in their operation and maintenance.

Operating Summary



Before proceeding please check the specification sheet at the front of this manual which options and accessories, if any, are fitted to your Priorclave. This will determine what you will need to read the instructions for these options later in this manual



1. Check electricity supply and that the power is switched to the wall socket
2. Press Door Hold button on the control panel. Hold it until 'Door Hold' is displayed in the temperature display. Wait for a short time until the temperature returns to normal, there is another bleep and the door indicator light can now be pressed again to release the lock.
3. Lift the locking handle until it reaches its safety stop. Do not lift handle until withdrawn as this may lead to damage.
4. Release the safety catch by pushing it forward with your thumb and lift the locking handle to the top of its travel. The door can now be opened.
5. Top up with water if necessary to just below the top of the weir.



ALWAYS CHECK THE WATER LEVEL BEFORE STARTING A CYCLE

6. Load the autoclave.
7. Set the temperature as required using the up/down keys.
8. Set the process time as required using the up/down keys.
9. Set / select other functions i.e. fast, rapid cooling etc., as required and if fitted.
10. Carefully close the door with the locking handle fully up.
11. Lower the locking handle to the bottom of its travel in one action to lock the door.
12. Wait a few seconds and start the cycle.

Cycle Abort and Thermal Lock Override

Aborting a cycle

To abort the cycle start by pressing the Door Hold button.

Thermal Lock Override

First abort the cycle as above.

After checking that there is no pressure within the autoclave turn the thermal lock key to the 'unlock' position.

Press Door Hold button keeping the thermal lock key held over.

Wait during the display until there is a bleep and the Door indicator light is on.

Keep the key held over until the display shows 'unlock'.


The key switch can now be released and the door opened as above.




If the key is released at any stage the procedure must be repeated to open the door.

Layout Diagram


Non Vacuum Models



Vacuum Models



Control Panel Details



- | | |
|---|---|
| 1 Door Button & Indicator | 7 Time Display & Setting Buttons |
| 2 Cycle Progress Display | 8 Program Buttons & Indicators (Option) |
| 3 Vent Button & Indicator | 9 Printer (Optional) |
| 4 Start Button & Indicator | 10 Setting Lock key switch (Optional) |
| 5 Option Setting Buttons & Indicators | 11 Thermal Lock key switch |
| 6 Temperature Display & Setting Buttons | 12 Hidden 'till lit f a |

Operation



Before using your Priorclave for the first time check that the supply is switched on that the water supply (if required) is available and that the Emergency Stop Button (if fitted) is not pushed in.



Before proceeding please check the specifications supplied with this manual to establish which options and accessories, if any, are fitted to your Priorclave. This will help you will need to read the instructions for these options later in this manual.

Opening the Door

1. Switch on the power at the wall socket. All indicators will light momentarily and a sounder will beep. This is to enable the indicators to be checked.



If you are powering up the autoclave for the first time, or if the autoclave is fitted with conductivity water detection and low water condition the error F00 will be showing in the temperature display. You may ignore this at this stage as opening the door and filling it with water will reset this error code.

2. The start and/or low water indicators should now be lit. Check that the pressure gauge is at zero and then press and release the door button, which will beep, and wait for a short time (about 20 seconds) until the door indicator illuminates and the sounder beeps again.
3. During the waiting time the temperature display will show 0.0 confirming that the autoclave is waiting during its safety delay.
4. On pressing the door button a second time you may hear the locking solenoid operate. You may now lift up the locking handle to the safety catch position (the door button will light when the door button is pressed and remain lit after the handle is lifted).



Do not attempt to open the door before the lock has released or damage to the locking mechanism may result which will not be covered by the warranty.

5. The door is retained by the safety catch. In the unlikely event of undetected residue being present the safety catch prevents the door from being blown open by breaking the seal whilst retaining the door, thus allowing the pressure to escape safely from a gap at the bottom of the door.
6. The safety catch is released by pushing it forwards with your thumb and then lifting the door handle to its fully open position. The door can now be pulled open gently.



Care should be taken when opening the door as it will be hot and steam may be present. Heatproof gloves and a face shield should always be worn when unloading an autoclave.

Checking Water Level.



ALWAYS CHECK THE WATER LEVEL BEFORE STARTING A CYCLE

The autoclave uses immersion heaters in a reservoir of water below the lower shelf to generate steam. The heaters are protected from boiling dry by a low water cut.

The reservoir should be filled with water to a level just below the top of the front of the autoclave.

If your autoclave is fitted with an automatic filling system then the autoclave will be filled with water to the correct water level as the door is opened.

Protection against the autoclave boiling dry and damaging the heating elements in low water in the autoclave vessel is provided by one of two methods:

Standard Protection

A temperature sensor attached to the heater in the autoclave chamber is connected to Tactrol controller. The Overheat temperature is factory set.

If an overtemperature is detected from the sensor then power to the heating circuit is switched off and fault code F018 will be shown in the temperature display.



Should this occur, then the autoclave should be switched off and allowed to cool temperature and pressure before opening for examination of the water level.

Conductivity Detection Protection

This system is ideally suited for use with softened water and will also be fitted if you fitted with an automatic filling system.

It uses a water conductivity probe (or 2 in the case of the filling system) which is using the conductivity of the water.

If the water level falls below the sensor the autoclave shuts down and warning is lit and fault code F004 is shown in the temperature display.



In hard water areas softened water must be used to prevent scale from forming. If manually filled then distilled water can be used but when the autoclave is new, may need to be added until the low water lamp is extinguished as the low water cut out operated by the water's conductivity level is always recommended before commencing a cycle.



Great care should be taken to ensure that the area of the low water sensor is clear as a build up of contamination here will prevent the low water cut out from working and do lead to heater damage.



The F018 or F004 indication must be manually reset when the door is opened by turning the program setting lock switch from position 1 to 3 and back again. If programmer is not fitted by pressing the Control Reset switch at the back of the autoclave. However the autoclave should be checked by an engineer for faults as soon as possible and cleaned or repaired if necessary before reuse.

A final level of temperature protection is provided by a thermocouple in line with the heater.

Loading

The autoclave can now be loaded with the items to be sterilised onto the shelves or baskets, or in the case of waste loads which may leak liquids when autoclaved into containers.



Care should be taken when loading the baskets or containers not to pack them too material. Ample room must be allowed for steam to penetrate the load properly for sterilisation will not be achieved. When using autoclave bags these should be left top of the bag rolled outwards, exposing the load to the steam inside the pressure.



Care should also be taken that the contents of bags and containers are not able to spill into the body of the autoclave vessel. Any such spill could block pipes and valves and will not be covered by the warranty.

For waste loads, which may leak liquids when autoclaved, watertight discard containers are strongly recommended.



Tests have shown that the depth of perforated discard containers should be no greater than 100mm. Prior to autoclave.

Settings.

Once the autoclave has been satisfactorily loaded the controls should be set for that you require.

If your autoclave has a setting lock fitted this must be set to allow the parameter to be altered.

Sterilising Temperature & Time Settings.



Research carried out by the Medical Research Council has recommended the following temperatures and times as being sufficient for complete sterilisation in autoclave

- 126°C for 10 minutes.
- 121°C for 15 minutes.
- 115°C for 30 minutes.

These temperatures and times relate of course to load temperatures and the aim cycle should be to achieve one of the above criteria in the coldest part of the load however are sensitive to elevated temperatures for prolonged periods, making full of the above impractical. However the disinfection of a full load, without necessarily reaching full Sterilising Temperature, is usually sufficient for most purposes.

Should you require a more precise method then the Load Sensed Process Timing be of assistance for certain load types. Autoclave is fitted with this option please refer description later in this manual.

Since there is a time and temperature 'lag' for a load, this should be compensated for either by increasing the process time, or by including in the cycle a period of free steaming with the vent closed to assist greatly with reducing the temperature lag.

This can be achieved by pressing the vent button manually and releasing it manually. Alternatively, the Automatic Free Steaming function will perform this function automatically if the option is selected.

If you have an interest in any of the options mentioned above, which can quite easily be fitted, please contact Priority Service America.

In conclusion, when setting up the autoclaving cycle a large safety margin should be within the settings.

More precise settings can be assessed by

Setting the process time

The process timer can be set to a time up to 999 minutes: the time required is set by simply using the time up/down buttons. Once the set time is displayed until the set temperature is reached, the process time begins counting down to zero in increments of one minute.



Setting process temperature

Pressing either the up or down button momentarily cause the current set temperature to be displayed. Subsequent use of the up/down buttons changes the set temperature. If no key is pressed for a short time, the display returns to showing the chamber temperature.



Selecting other functions

The function select keys may be used to switch the Free Steaming, Cooling, Media Warming Option, and optional functions such as Load Sensing, Auto Vacuum cycle on or off at any time other than when a cycle is running. An indicator illuminates when a function has been selected. If an option is not fitted (or permitted in the selected program memory models) pressing the appropriate key will result in a visual and a fault being signalled and the function will not be selected.

Automatic Timed Free-Steaming

What is free steaming?



Incorporating a period of free steaming into a cycle can improve air removal from difficult loads and/or reduce temperature lag between the load and the autoclave, reducing process times and higher temperatures. Free steaming introduces a stage during heating up to Process Temperature, when a solenoid valve at the rear of the autoclave is opened for a set time. The valve opens at a factory set temperature of just above 100°C and is held open for the time detailed below. During this time steam is being generated in the chamber in large quantities, this creates turbulence that passes through the load before escaping through the valve, the turbulence that can assist with air removal.

Setting the free steam time.



If free steaming is required this is selected by pressing the Free Steam button. The indicator lights up to show that Free Steaming is selected. The time display will now flash indicating the Free Steam time, not the process time, is currently being displayed. The time can now be set (in minutes) using the up/down buttons. If no further changes are required a short time the display stops flashing, and reverts to showing process time. If you wish to change the Free Steam time or make further changes the Free Steaming should be deselected, then reselected.



Caution should be used before setting the Free Steam time longer than 15 minutes. Excessive free steaming times can use a large amount of water, increasing the possibility of completing due to a Low Water condition.

Free Steam temperature setting

Timed free steaming will commence at a temperature slightly above 100°C, which has been set at the time of manufacture. If required, this temperature may be increased by qualification and the turbulence caused by the escaping steam pressure can further assist with air removal, it is desirable however to connect the autoclave to a drain and vent pipe in accordance with the Installation manual for this model to avoid the release of pressurised steam into the laboratory.

Performance can be improved even further by fitting the optional pulsing system.

Pulsed Free-Steaming



(Optional Function included for autoclaves with vacuum cycles)

With certain loads and in certain situations the efficiency of the free steaming process can be improved by pulsing. When available according to program free steaming commences as described above but at a higher temperature. Instead of remaining open for the free steaming period the vent valve shuts off at a lower temperature. The autoclave returns again to the temperature at which the valve opens again. The autoclave will continue this cycle for the time set when selected. This continual pulsing of steam out of the autoclave creates considerable turbulence within the autoclave, helping to draw air from the load.

If fitted, this function is program specific and is usually set up during commissioning on particular programs in response to customer requests at the time of ordering.

When selected as part of the program the pulsing function is standard steam function described above.



Setting of the free steam time for a particular program is as described above.

Pulsed free steaming is not suitable for bottled liquids and should not be selected for these type of load.

Load Sensed Process Timing



(Optional Function)

Load Sensed Process Timing Function

If this option is fitted, the autoclave will be provided with an additional probe inside the autoclave chamber. This probe can be positioned in the load, ideally in the coolest part of the load. When this option is selected, the autoclave will heat to the set chamber temperature. However, when the set temperature is reached the process time will not begin to count until the load temperature as sensed by the additional probe reaches a temperature just below the set chamber temperature. The cycle will then proceed in normal manner.

Load Sensed Process Timing Purpose

The use of load sensed process timing can greatly assist the sterilisation of certain types of difficult dense loads, such as large baskets of bottle caps, pipette tips or animal material that the load reaches set temperature. The system is also very effective for bag loads, however as these tend to melt down around the probe, consumption of probe is high. For this type of use load validation may prove to be more successful and effective in the long term.

Positioning Load Sense Probe

The probe should be positioned in what is anticipated to be the slowest part of the load to heat up, for example the centre of a large densely packed load, or the largest of a group of bottles. This is important, as there may be large variations in temperature distribution throughout the load.

Temperature variations can be reduced by the use of timed free steaming.

Load Sensed Temperature Probes

The load sensed process timing option utilises a temperature probe connected directly to the main processor board.

Replacement temperature probes are available from Priorclave.

Rapid Cooling



A fan is fitted into the bodywork of the autoclave to direct cool air over the autoclave chamber.

If selected by using the cooling button, the cooling fan will switch on automatically during the cooling stage of the cycle. There are three possible settings for rapid cooling, as follows:

Off- No indicators lit.

Immediate start The cooling fan does not operate at all during the cycle.

Left hand indicator-lit press of the cooling button.

The cooling fan starts as soon as the cooling stage is reached.

Delayed start Both indicators lit press of the cooling button.

The cooling fan starts after the autoclave chamber has cooled. This setting is useful when autoclaving some fluid loads, as bringing the chamber on at temperatures above 100°C may reduce the chamber pressure too rapidly, causing the load to boil.

In both cases the fan will switch off automatically when the cycle has completed the stage.

Media Warming



If this highly useful feature is selected the autoclave will cool to a temperature of 45°C. The temperature will then cycle between approximately 45°C until the door is opened. This allows, for example, nutrient media to be held as a liquid until it is needed especially when used along with the delayed start function

Delayed Start Time

The autoclave can be set before a cycle to start at a time

To access these settings turn & hold the thermal lock key in the override position. Press up or down keys. Release the thermal lock key, the time displayed on the temperature display, 10 will be displayed on the time display. The temperature display now shows the list of operating parameters, the value for the parameter is shown in the time display through the list of available parameters using the temperature up/down keys.

After no keys are pressed for eight seconds the display returns to normal.

The function of these settings is as follows:

Temp. Display	Time Display	Function	Action
1	0-24	Delayed Start Time Hour	Enter required Start time hour (24 hours)
2	0-59	Delayed Start time Minute	Enter required Start time minute
3	0/1	Start Delay Set On/Off	0= OFF 1= ON
+ The time is set in real time, therefore the clock has to be correctly set for this properly. After one delayed start operation, delayed start automatically switches off, and returns to normal operation			

For instructions on setting the clock time and for other operator settings please see section Changing Date & Time in this manual.

Vacuum Options (MVA (vacuum) models only)

(Optional item)



It is strongly recommended that to obtain optimum performance from Priorclaves fitted with vacuum options that commissioning and/or load validation tests are carried out by a Priorclave engineer. If no particular programs have been specified your autoclave will be set with the following programs:

- Program 1: Non vacuum Cycle (Pre-Cycle Vacuum can be selected)
- Program 2: Pre-Cycle Vacuum and Vacuum Cooling
- Program 3: Pre-Cycle Vacuum and Vacuum Cooling
- Program 4: Non vacuum Cycle (Pre-Cycle Vacuum can be selected)
- Program 5: Non vacuum Cycle (Pre-Cycle Vacuum can be selected)



PreCycle Vacuum

The precycle vacuum is selected using the function select key on the control panel and hand indicator lit the Cycle Vacuum is selected. With the Pre Vacuum selected a vacuum pump will run at the beginning of the cycle, removing much of the air from the load and chamber. At a preset level of vacuum the control system switches off the pump and the cycle begins. By default two vacuum stages will be performed with a heating stage in between. Pre cycle vacuum is essential when autoclaving loads containing densely packed porous materials.



Vacuum Cooling Suitable for Non-Liquid Loads Only

A vacuum cooling cycle can be selected by means of the function select key. With the option selected the right hand lamp will illuminate. When this option is fitted it can be run separately from a Cycle Vacuum. With the option selected, at the end of the process time the autoclave vent is opened and the autoclave cools to a set temperature with the cooling fan(s) operating. When the set temperature is reached the cooling fan(s) stop and a partial vacuum is drawn. This has the effect of evaporating liquid on the load to cool rapidly. After a preset time air is admitted to the vessel and this process is repeated a number of times. At the end of this stage the autoclave passes immediately to cycle complete.



Post cycle vacuum cooling must not be selected if the load contains bottled liquids, regardless of how these are contained. All liquids in the load will be contained, contaminating the inside of the autoclave chamber. Sealed containers of liquid are likely to explode. Unexploded containers may be in a dangerously unstable condition when removed.



Drying Cycle Suitable for Non-Liquid Loads Only (Optional Vacuum Options Fitted)



This option must be selected for attachment to a particular program in the control panel during commissioning.

A drying cycle can be selected by means of the function select key within a program designated as a drying program. With the option selected the right hand lamp will illuminate. When this option is fitted it can be run along with or separately from a Cycle Vacuum. With the option selected, at the end of the process dwell time the water charge is drawn from the autoclave, and the autoclave cools to a set temperature. When this temperature is reached a partial vacuum is drawn and heaters on the outside of the autoclave vessel are switched on. This has the effect of evaporating liquid on the load. After a preset time air is admitted to the vessel and this process is repeated a number of times. At the end of this stage the autoclave passes immediately to cycle complete.



Post cycle drying must not be selected if the load contains bottled liquids, regardless of how these are contained. All liquids in the load will be contained, contaminating the inside of the autoclave chamber. Sealed containers of liquid are likely to explode. Unexploded containers may be in a dangerously unstable condition when removed.

Multi Program Memory Options

When this option is fitted, five program number keys are provided to the right hand panel, each with two indicators. The indicators on the left are for programs 1 to 5 and on the right for programs 6 to 10. If the Priorclave has been specified with a 10 program memory only the first five programs will be active. As each program is selected the indicator illuminates and the previously selected indicator is cancelled. Pressing the button toggles between the two program numbers shown on the button.



When the program memory option is fitted, position setting lock switch is fitted. These setting positions allow different levels of access to settings as follows.

Position 1. Only the currently selected program can be run.

Program settings cannot be changed.

Position 2 All programs can be selected and run.

Program settings cannot be changed.

Position 3. All programs can be selected and run.

Program settings can be changed freely.

NOTE: The setting lock key can only be removed in positions 1 and 2

Programming of settings is the same as with the machine, but the required program number should be selected before setting. The settings entered can then be re-used by simply reselecting that program number.

Closing the pressure door


When you have set up the cycle, close the door with the locking handle in the raised position. Then, in a single action push the handle down into the locked position. The door will now be properly secured.

Starting a cycle

Ensure the door is properly secured and the indicator is illuminated. To start the cycle, press the start button.

The first segment of the cycle status indicator bar will illuminate and the autoclave will gradually heat up to process temperature. The cycle status indicator will also go through its stages to give 'at a glance' information on the cycle progress.


Once a cycle has been started the function selection settings cannot be changed. To do so will cause a fault to be signalled. If changes are required, the cycle must be stopped by pressing the start button again.

Pressing the Free Steam Button  during the Free Steam part of the cycle will give an indication of Free Steam time remaining as a flashing display in the timer window.

On Priorclaves not fitted with the setting lock key, changes can be made to the program time and temperature settings once a cycle has started. At the end of the cycle, the settings will reset to their original settings.



When the autoclave has reached its operating temperature, the door cover and other near surfaces may become uncomfortably hot. Take care touching these surfaces without protection.



Vent button.

The vent button may be used at any stage during the cycle. When used it opens a valve at the back of the machine. It may be left open for free steaming to achieve better penetration of the load. Automatic Freesteaming has not been selected. Care should be taken with this manual method however as failing to switch off the vent will waste a large amount of steam and to eventual cycle failure. It may also be used with certain loads as a means of venting the autoclave. If used for manual steaming the vent button must be manually released before pressure will build up and process temperature can be achieved.



Care should be taken if using the vent button when the autoclave is pressurised. autoclave under these circumstances with a liquid load may lead to the load boiling and glassware may be broken.

During the process time.

Once set temperature is reached, the process time will begin to count down and the indicator will illuminate. If the Load Sensed Process Timing Option is fitted and selected, there may be a delay between the autoclave having set temperature and commencement of the process time whilst the load reaches set temperature.

During the process time a check should be made that there is correct correlation between temperature and pressure readings on the control panel. A steam indicator is located at the back of the autoclave. The check should be made to ensure that air has been purged from the autoclave. Generally, a pressure reading higher than would be expected to indicate entrapped air in the autoclave.

If for any reason the temperature is forced outside range or power to the autoclave is removed during the process time, the cycle will abort and the fault indicator will show a fault code of either F005 or F006. This is to ensure that loads that have not been subjected to the required cycle parameters are not processed incorrectly. The fault condition is cancelled by:

If no setting lock is fitted:

pressing the reset button on the top panel on the right hand side of the control panel

or

If a setting lock is fitted:

turning the setting lock key to the enable position and then to the disable position. If the autoclave was in the enable condition when the fault occurred, then it must first be turned to the disable position.

Power loss during a running cycle

Loss of power to the autoclave whilst the cycle is running can have one of two consequences depending on the stage of the cycle at which the power loss occurs.

If the power loss occurs at any point before completion of the sterilising stage of the cycle (process time) then the error code E004 will be shown in the temperature indicator. The cycle will be stopped. This is to ensure that loads that have not been subjected to the required cycle parameters are not assumed to have been processed correctly. The fault condition is described above.

If the power loss occurs after completion of the sterilising stage of the cycle (process time) the cycle will stop but error code E006 is not shown in the temperature indicator. The autoclave was cooling or running a post vacuum process at the point of power loss. The cycle will resume after the power loss.

In many cases and with care this distinction between the stage at which power loss occurs can be helpful. It allows operators to determine whether the load is safe to be disposed of rather than having to be run through the autoclave again.

Cooling

After completion of the process time the autoclave moves into the cooling part of its cycle: this is shown on the cycle status indicator. Cooling has been selected this will be switched on automatically according to the cooling strategy selected. It will use convection.

If there is a power fail while the autoclave is in the cooling part of the process it will resume in the cooling phase once power is restored.

If the power should fail before sterilisation is complete then the cycle will be stopped.

Thermal lock

Under normal circumstances the autoclave cannot be opened until the temperature simulator probe, which has a cooling rate assimilated to a bottle of fluid, has fallen to a point where the yellow thermal lock status indicator will illuminate. The temperature by the temperature indicator will be significantly below 100°C as the probe measures the temperature of the open chamber space. Pressing the door button before the thermal lock has a fault to be signalled. The temperature at which the thermal lock operates is factory set but must only be done following commissioning by qualified personnel. It can be overridden using the key switch on the control panel. The keys for this switch are provided in this manual.



Overriding the thermal lock will cause the main vent to open. Great care should be taken when using the key switch since liquid loads could boil over if vented at elevated pressure and glassware could be damaged.

There are circumstances, however, when quicker access to the load is required. When necessary, first abort the cycle by pressing the start button. Then turn the key in the key switch to the 'on' position and holding it in this position press the door button and wait while the thermal lock message is displayed until the door lamp illuminates. Finally press the door button to release the door lock. The thermal lock key can now be released. If the key is released before this message is displayed the display will not reset and the autoclave cannot be opened. To reset the display, repeat the above procedure and open the autoclave.



Great care should be exercised when using the Thermal Lock Override, especially with liquid loads. Even at temperatures below 100°C a liquid load in sealable glass containers will not be safe. For the above reasons responsible personnel should keep the Thermal Lock Override in a safe place away from the autoclave to prevent access to it by unauthorised personnel.



Under certain cycle abort or failure conditions the thermal safety lock can latch in condition. This is because the control system will always go to the safest condition any uncertainty about the cycle end circumstances. To override simply go through the door open or close procedure using the thermal lock override key. Operation will return normal as soon as the next cycle is completed satisfactorily.

Cycle complete

When cooling to the 'thermo' complete, the thermo indicator will illuminate, and the autoclave will emit a bleep for a short time (about 10 seconds). If the Cooling System is selected it will automatically switch off at this point. The autoclave door will then open and unload.

Media Warming

If this has been selected, the autoclave will remain at the temperature after the cycle complete, until the door is opened or the cycle otherwise aborted.

Opening the autoclave to unload and for the next cycle is simply a repetition of steps 2.

Aborting a Cycle

On occasions it may be necessary to abort a cycle before its completion. In order to do this simply press the start button.

Emergency «E» Stop Button (if fitted)

In an emergency, pushing in the emergency stop button will cut the electricity supply to the heaters and control system.

After use this must be released using the key provided.

This key should be kept by responsible personnel in a safe place away from the autoclave. Access to it by unauthorized personnel should be prevented.

Operation with Options & Accessories

The following descriptions detail how to operate and gain maximum benefit from accessories that may be fitted to your Priorclave.

Setting LockKey-switchOption

Fitted on Priorclaves without program memory to give an optional level of security. It has two settings only, which are equivalent to positions 1 & 3 described previously. The key can only be removed in position 1.

Printer

The printer if fitted is mounted on the right hand side of the control panel. This provides a record of the cycle as well as an indication if any faults have occurred. The information is printed as follows:

USER NAME (if provided at time of ordering)
DEPARTMENT (if provided at time of ordering)
AUTOCLAVE SERIAL NUMBER
DATE (in the format dd/mm/yy)
CYCLE NUMBER
PROGRAM NUMBER (if multi-program memory option fitted)
TEMPERATURE AND TIME AT CYCLE START (time initially set to G.M.T)
TEMPERATURE AND TIME AT END OF FREESTEAMING
TEMPERATURE AND TIME AT START OF PROCESS TIME

The temperature and time are then recorded at time intervals, until the end of the process time.

TEMPERATURE AND TIME AT END OF PROCESS TIME
TEMPERATURE AND TIME AT CYCLE COMPLETE.
CYCLE PASS/FAIL/ABORT


(Fail means that a fault signal has occurred during the cycle or that the cycle has been aborted. See-Warning Indicators)

The above information will be printed in the order listed, allowing the information to be read as it is printed.

Power On Self Test

The self test procedure is initiated by applying power to the printer while the paper feed button is depressed. When the paper feed button is released a test print will be produced.

Replacing Paper Roll



Correct Paper Path


Access to the printer for changing the paper roll is provided by a hinged front panel. The catch on the left side of the printer front panel should be depressed to release the front panel.

The paper roll carrier is mounted on the rear of the printer and will swing out from the main body for ease of roll changing.



Always use the correct paper. Suitable paper is available from Priorclave Service or your distributor.

Changing Printer Ribbon



Ribbon fitting detail

With the door of the printer open, push the printer door at the point shown in the diagram. The outer door cover should come away exposing the print ribbon mounting plate above.

The old print ribbon can be lifted out and a new one put in its place. If necessary, slack on the ribbon with the tensioning wheel.

The paper should pass between the ribbon and the bottom of the ribbon cassette. The two parts of the door can now be pushed together and will click into position.

Serial Interface

The autoclave is supplied with an externally mounted serial connector and a suitable cable for serial or USB connection to a computer. Also supplied is a website address from which to download the latest software and a full operation manual for the system.

USB Interface

The autoclave is supplied with an externally mounted USB port. Also supplied is a website address from which to download the latest software and a full operation manual for the system.

Automatic Waterfill Option

(Not suitable for RO, or Purified water supplies)

This option consists of a water tank mounted on the side extension of the autoclave. It contains a conductivity water level probe set to detect minimum water level and a second conductivity water level probe set to detect water level in the fill level for the autoclave chamber.

The water tank is fitted with a float valve to automatically control fill level.

Whilst the autoclave door is open during loading and unloading a valve is opened and water flows into the autoclave vessel from the water tank until the upper water level probe is triggered.

This system is additionally fitted with a timer that automatically shuts off the water supply after a preset time. This is fitted as a measure to prevent flooding in case of a water level probe failure.

If the timer times out before the water level reaches the probe, a 'Water Fill Fault' is displayed on the control panel.

To start filling again, reset the fault (see Warning Indicators and Fault Codes). The probe should be regularly checked and cleaned (see Maintenance Section for details).

During normal operation only small amounts of water will need to be replaced, however, first filling the autoclave or whilst filling after cleaning, more time will be required.

Air Intake Filter

(Optional or fitted as part of a vacuum drying system)

When fitted to the autoclave this system ensures that air drawn into the autoclave during the cooling stage of the cycle is first passed through a bacteriological air filter. This is located at the back of the autoclave.



The filter should be regularly changed to maintain its effectiveness

Vent Filter Option

When fitted to the autoclave this system passes all autoclave discharge through a vent filter housed inside a pressurisable stainless steel housing. At the end of the cycle any unfiltered condensate from the filter housing is returned to the autoclave. The correct operation and effectiveness of the filter system should be regularly checked as described in this manual.



Internal Validation System

If fitted this system continuously monitors the performance of the temperature control system with reference to an internal reference standard. If any problem is detected the autoclave is stopped and a fault is signalled. (See Warning Indicators Codes for details.)

Accelerated media cooling

Unless specified otherwise before delivery of the autoclave or set up during the commissioning process the Media cooling function is set to operate by default program 1.

It can easily be set to operate on different programs by a Priorclave Engineer or your approved Priorclave representative.

When selected the autoclave will carry out a standard media cooling program as normal with the cooling fan coming on at the end of the process time.

As the autoclave cools below 100°C the pressure on the gauge will gradually drop below atmospheric pressure. This assists greatly in the cooling of liquid loads giving up to 50% of the total cooling time to cycle completion at a safe thermal lock temperature.

Once the cycle has completed any vacuum in the autoclave is automatically released and the door can be opened in the normal way.

For best results larger sized bottles should be used, i.e. 500ml volumes in a 1000ml bottle.



Delayed start cooling (both lights lit on the cooling button) may be required if excessive volume losses during cooling are observed.

Chart Recorder/ Data Recorder / Data Logger

If fitted, the recorder power input is connected directly to the autoclave in such a way that the recorder will only operate during the autoclave cycle, i.e. from the pressing of the cycle start button to the cycle complete. At this time power to the recorder is cut off. Unless otherwise specified, single channel units record the temperature of the temperature probe and in the case of two channel units the second channel records the temperature of the load probe.

For more details on individual recorder function and operation please refer to the recorder or logger manufacturer's manual supplied with the recorder.

Changing Date & Time

A number of additional control system settings are available. To access these settings hold the thermal lock key in the override position. Press the up or down keys. Release the thermal lock key. The value displayed on the temperature display, default 0 will be displayed on the time display. The temperature display now shows the list of operating parameters, the value for the parameter is shown in the time display through the list of available parameters using the temperature up/down keys.

After no keys are pressed for eight seconds the display returns to normal.

The function of these settings is as follows:

Temp. Display	Time Display	Function	Action	
		The autoclave can be set for the cycle to start-programmed delay, for example to allow a media preparation cycle to complete shortly prior to the working day. Setting the value of parameter 1 to 1 in the time display will delayed start on.		
1	0-24	Delayed Start Time Hour	Enter required Start time hour (24 hr clock)	
2	0-59	Delayed Start time Minute	Enter required Start time minute	
3	0/1	Start Delay Select On/Off	1= ON 0= OFF	
*	4	0-999	Print Interval	Enter time (minutes) between printing process time (0= printer disabled)
5		Year Setting	Enter Year	
6		Month Setting	Enter Month	
7		Date Setting	Enter Day of month	
8		Hour Setting	Enter Hour (24 Hr Clock)	
9		Minute Setting	Enter Minute	
10		Second Setting	Enter Second.	
		Scroll back up to parameter 1 to confirm the new or current time		
#	11	1-999	Cycle Repeats	Enter Number of Cycles Required
+	The time is set in real time, therefore the clock has to be correctly set for this properly. After one delayed start operation, delayed start automatically switches off, and returns to normal operation.			
#	Models fitted with optional Cycle Repeat Facility only			
*	Models fitted with 5 or 10 Program Memory Setting marked are program number related, and therefore should you wish to use different values for these in different programs this can be done by changing the program number when the correct program is selected.			

Warning Indicators and Fault Codes

If one of these faults occurs please contact your service provider with details of the serial number of your autoclave



Do not attempt to rectify these faults (with the exception of the Service indicator and the water fault) yourself.

On the control panel there are a series of indicators which will appear in conjunction with a fault code in the test panel. The meaning of these warnings, why they appear, and what to do when they appear, is as follows

SERVICE

This means that 500 cycles, or six months have passed since the autoclave was last serviced. An engineer will cancel the message when the autoclave is serviced.

WATER + FAULT CODE F004

The water level has fallen below the minimum level and must be topped up before the autoclave can be run. The warning will automatically cancel when the door is opened and the water is topped up. The low water condition may have caused a running cycle to abort, and the autoclave will need to be autoclaved again.

O/HEAT + FAULT CODE F003

If fitted, the heater over temperature protection temperature sensor may have sensed that the heating element became too hot. This is probably due to a low water condition, which is sensed by the low water probe. The water level and the condition of the probe (see section 4.1) should be checked before attempting to use the autoclave again.

If heater over protection is not fitted then the heater cut out will only operate under exceptional conditions, such as a failure of the temperature control system. The next attempt to run the autoclave should be closely observed and if problems persist contact Priorclave Service.

FAULT + FAULT CODES F000, F002, F005, F006, F007, F008, F009, F010 & F011

The fault indicator illuminates under conditions that may invalidate the autoclaving process. The autoclave may result in the load requiring to be autoclaved again. The fault code is triggered by any of the following:

- F006 Power to the autoclave being interrupted when a cycle is in the heating or processing stage of the cycle.
- F005 The chamber temperature falling below the set temperature by more than the process dwelltime.
- F002 Failure of the temperature control, display, or load simulator temperature sensor.
- F000 If your autoclave is fitted with the optional load simulator system, an error in the temperature measurement system is signalled by fault code F000. Usually this would mean that a critical fault has developed in the temperature measurement system, however, as the detection system is extremely sensitive it is possible that it may be triggered by fluctuations in the electrical supply. If fault code F000 appears it may be cleared by the method described below. If the fault code will not clear, or continues to appear then the user cannot correct the fault. In such cases please contact Priorclave service or your local Priorclave approved service agent.
- F007 Vacuum stage timeout (loop break). The autoclave has not achieved the set vacuum level during the process cycle vacuum stage during the process time.
- F008 Heating stage timeout. The autoclave has not reached process temperature within the set time.

- F009 Vacuum cooling point not achieved. The autoclave has not achieved a low enough vacuum during the post cycle vacuum stage (Vacuum Cooling or Drying Cycle)
- F010 Air detector input activated. If fitted the air detector system has over pressure condition symptomatic of excess air remaining in the load.
- F011 Printer Timeout / Malfunction. The control system has not received confirmation from printer within its preset timeout.
- F012 Door microswitch fault. If a door microswitch opens during a cycle this fault code is displayed.
- F013 Jacket Timeout If a jacket is fitted it has not reached the required temperature within the set time. This would indicate a problem with steam supply or inlet or drain valve operation.
- F014 Jacket Over temperature If a jacket is fitted the temperature has exceeded the set temperature.
- F015 Jacket under temperature If a jacket is fitted the temperature has fallen below the pre-set operating temperature band.
- F016 Water Fill Timeout The upper level water probe level has not been reached within the set time for filling and the filling operation has been stopped. This function prevents continuous unsupervised operation of the water fill, which could lead to flooding.
- F017 FreeSteam-During Pulsed Free Steaming operation the lower of the two set temperatures has not been achieved. The temperature has not fallen sufficiently following the opening of the steam valve.
- F018 Heater Overheat. If this fault occurs the most likely cause is a Low Water condition. Check the water supply is turned on and the condition of the heater before resetting this fault.

LOCK

This warning will light when the thermal lock is in the override position.

LOAD + FAULT CODE

This warning is activated in the event of the failure of the load sensing temperature sensor. The temperature sensor should be replaced as soon as possible. Great care should be taken that loads which would ordinarily be autoclaved with dry process timing are adequately sterilised.

CANCELLING FAULT MESSAGES

The fault messages are cancelled by first correcting the source of the original fault and then setting the lock key switch to position 3. If a key is not fitted they are cancelled by pressing the reset button.

If 2 or more faults occur at the same time, the one with the highest priority is displayed first (the highest priority and F012 is the lowest.) If a higher priority fault is cleared it will display the next active fault, unless this too is cleared by the same action.

Routine Operator Maintenance



Before carrying out any maintenance work check the autoclave for any visual signs of contamination or damage. Should any contamination be apparent contact the relevant person of authority before proceeding.



Do not attempt to carry out any work unless you are competent to do so.



Disconnect or isolate the machine from mains power supply before removing any parts commencing any maintenance work.

Ensure that any electrically locked doors are open before disconnecting power.



In the event of any difficulty or doubt about any maintenance or service procedure contact Priorclave Limited or your nearest Priorclave approved service agent immediately.



For full maintenance and servicing details please refer to the Maintenance Manual.

Daily Maintenance

LOW and FILL Level Water Level Probes (where fitted)


To ensure protection from boiling dry, the insulated low water probe between the stainless steel tip and pressure vessel wall should be scrubbed clean to prevent from being short circuited (see diagram). The sensor also be kept clean to ensure good contact.

Water Level and Condition



The standing water within the vessel should be regularly removed, ideally with a wet/dry type vacuum cleaner to prevent the build up of spilled media and potentially corrosive chemicals. With the vessel emptied of water, the heating element(s) should be wiped with a damp cloth to remove build-up of Limescale.

If destruct type loads are regularly being processed consideration should be given to the use of suitable containers for such loads to minimise spillage. A range of such containers suitable for individual Priorclave models is available from Priorclave.



Weekly Maintenance

(To be carried out in addition to your daily maintenance program.)

Autoclave Cleanliness

Check exterior of machine and the inside walls of the pressure vessel for cleanliness, particularly around operating parts and external switches and dials. Use antibacterial wipes to clean exterior panelling.




Under no circumstances should an abrasive or chemically aggressive cleaner be used on the pressure vessel. Use of chlorine or hydroxide based cleaners is not recommended (see notices at the beginning of this manual).

Door Gasket (s)

To prolong the life of the sealing gaskets it is advisable to lubricate the sealing face, carried out with the pressure door in the open position by applying high melting point grease to the exposed surfaces of the gasket, after cleaning the gasket and inspecting for damage. High vacuum grease is ideal for this purpose. Ensure that the grease is spread evenly around the surface of the gasket.

Door Closure Arms

To ensure a free action of the Closure Arms, they should be lubricated and free from dirt. The mating face of the Arm should be lightly coated with Copper Grease to ensure smooth action.



Drainage

If the autoclave has been plumbed directly into a drain using flexible tubing, this should be checked for any signs of blockage, obstruction or damage. Also ensure that both tubes are connected as originally intended. Check for any obstruction to the safety valve which must remain exposed and unconnected to any form of drainage at all time.

Indicator Lamp Check

To ensure that the Control Panel is displaying correct information, perform a check of the function of the Control Panel by tripping the Autoclave power supply. The process of starting will illuminate ALL Control Panel lamps for a period of approximately 3 seconds. If any lamps are seen to be OFF during this time, report the matter for investigation.

Monthly Maintenance

(To be carried out in addition to weekly maintenance program.)

General Cleaning

Check exterior of machine and inside pressure vessel for general cleanliness, particularly operating parts and external switches and pins. Wipe over all surfaces using antiseptic wipes.

Empty all water reservoirs using a vacuum water remover if available. In order to prolong the lifespan of the heating elements, water level probe, etc., it is advisable to avoid the build up of media, chalk, limescale, etc. on these parts or in the reservoir area. This build up is minimised by the use of softened water and a discard load container.



Under no circumstances should an abrasive cleaner be used on the pressure vessel. The use of chlorine or hydroxide based cleaners is not recommended (see notices at the beginning of the manual).

General Operation

The general operation and performance of the autoclave should be observed frequently to ascertain consistency. Any fault or defect reported or rectified immediately, should be entered into the Maintenance Log. (This will assist the service engineer in locating any persistent faults and reporting it to the manufacturer.)

Automatic Waterfill

The condition of the float valve tank and its various connections should be periodically checked and cleaned if necessary. Particular attention should be paid to ensuring that there are no restrictions to the tank overflow. In hard water areas the function of the float valve is of particular importance.

Routine Operator Maintenance

fitted to the water supply should be tested and the heating elements checked for scaling.

Bi-Annual Maintenance

(To be carried out in addition to weekly and monthly maintenance programme.)

Hinges

With the pressure lid in the open position the hinge should be cleaned and lubricated with melting point grease.

Checking Temperature Control and Pressure Gauge

During the DWELL stage of running autoclave cycle the Process Time is run for at least five minutes check the reading shown by the temperature display against that of the Steam Table manual.

A consistent disparity of readings would indicate either trapped air or a calibration issue. A note made of any observations in the Maintenance Log.



The pressure gauge and control system fitted to Priorclave are extremely reliable instruments such are unlikely to produce false readings. Therefore it is more likely that any deviation values given in the steam table is caused by incorrect purging etc.

Fault Finding & Rectification Guide



If one of these faults occurs please contact your service provider with details of serialnumber of your autoclave

Do not attempt to rectify these faults (with the exception the low water fault) y

Symptom	Possible Cause	Possible Solution
No Power	Power switched off at isolat	Check
	Circuit Breaker Tripped	Call Engineer
	Electrical Failure	Call Engineer
	Emergency Stop Button Pushed If (if fitted)	Release
Cycle does not commence when start button is pressed (Fault indicator flashes)	Door is not closed properly	Open & Reclose.
	Microswitch Failure	Call Engineer
Heating low or not apparent	Controller incorrectly set	Check setting & reset if necessary
	Circuit Breaker Tripped out	Call Engineer
	Heater(s) Failed	Call Engineer
	Heater(s) Over scaled	Call Engineer
	Fault in control circuit.	Call Engineer
Low Water Indicator (F004) stays lit when correct water level is achieved	Low Water Conductivity Distilled water	Add Tap Water
	Probe dirty, damaged or missing	Visually check probe. Clean necessary
	Wiring connections loose or damaged	Call Engineer
Autoclave does not pressurise	Vent button in open position	Check Vent Indicator and switch off
	FreeSteamtime not completed (if option is selected)	Check Cycle Progress Indicator
	Air purge valve failure	Call Engineer
	Vent Valve stuck open	Call Engineer
	Safety Valve stuck open	Call Engineer
	Door incorrectly closed	Check door position
Incorrect Temperature/Pressure Correlation	Air not fully purged from autoclave. Due to :	
	Incorrect Load Packing	Re-Load and restart cycle with (m FreeSteaming).
	Faulty Air Purge Valve	Call Engineer
	Water covering Controller	Check Water Level
	Faulty Controller or Gauge	Check function and calibration

Fault Finding

Symptom	Possible Cause	Possible Solution
Safety Valve opening	Temperature set too high	Check Temperature Setting
	Contactor Failure	Call Engineer
	Safety Valve Faulty	Check the seal on the valve is intact and has not been tampered with that the lifting handle is not obstructed. Check the pressure on the gauge is above 2.4 Bar. If not, call Engineer
	Output board Failure	Call Engineer
Door will not open once Autoclave Thermal Lock Temperature (WAIT) has cooled to 80°C	Thermal Lock is set to Load Chamber temperature yet reached in Load Simulation	not Chamber temperature.
Door does not open when door button is pressed	Safety delay not completed	Wait while = is displayed in the temperature display
	Thermal Lock Temperature (WAIT) yet reached in load simulation	Check Cycle Complete indicator in cycle progress display.
	Faulty Door Solenoid	Call Engineer
	Thermal Lock previously overridden	Use thermal lock key to open door. Normal function will return after complete cycle has been run.
Fault Indication will not go out	Fault not Reset. (Refer to manual for list of Fault Codes)	Call Engineer



Autoclaves are pressure equipment and as such are potentially extremely hazardous. They must be correctly serviced.

If you have any doubts or If you do not feel competent to carry out any of the above then do not hesitate to contact Priorclave North America on (800) 743-1590 or your nearest Priorclave approved service agent.

This device may be fitted with a resettable thermal cut out to protect the heater. In the event of failure of this device contact Priorclave Service or your local service agent.

Steam Table

Temperature		Pressure (BarG)
(°C)	(°F)	
100	212	0.00
105	221	0.20
110	230	0.43
115	239	0.69
120	248	0.99
121	250	1.06
122	252	1.13
124	255	1.25
126	259	1.35
128	262	1.55
130	266	1.70
132	270	1.86
134	273	2.04
136	277	2.21
138	281	2.40
140	284	2.60

Correct Correlation between Temperature and Pressure shows correct c
the autoclave and that air purging is satisfactory.



This table is accurate at sea level and at moderate altitudes. It will be suff
accurate for its intended purpose. However, at higher altitudes the pres
indicated will be slightly higher than those shown above.

Other Items Fitted

