













Illustrated Spares Price List

















Priorclave

13/01/2023

Part Number	Description	Price	
AFL/ACS/EXH	Containment Filter Element - External <i>For Filter Housing</i>	£471.25	
AFL/ACS/EXI	Containment Filter Element - Internal	£601.24	
AFL/ACS/SML	Filter - Air Intake - SMALL	£141.41	
CAS/MID/000	Castor Wheel Set - Compact C60 model <i>Set of 4</i>	£30.15	
CON/230/24A	Contactor - 50amp	£101.53	
CON/230/74A	Contactor-74 Amp 230V	£274.02	
CON/240/09A	Contactor - 9 Amp 240V	£49.87	
CON/240/23A	Contactor - 23 Amp 240V	£85.72	
CON/240/45A	Contactor - 100 Amp 240V	£189.92	
CON/240/811	Overload Relay (8A-11A Range)	£58.94	
CON/24V/09A	Contactor - 9 Amp 24V	£62.51	
CON/OLD/004	Overload Relay (2.4 - 4.0 Amp range) <i>Omron Type</i>	£44.11	
CON/OLD/69A	Overload Relay (6-9 Amp range) <i>Omron Type</i>	£78.39	

Part Number	Description	Price	
CON/OLD/OL6	Overload Relay (4-6 Amp range) <i>Omron Type</i>	£40.74	
CTR/ACS/001	Cycle Counter Panel Mount 12V DC	£32.69	
DTM/240/002	Delay Timer - 2 Contact	£368.35	
FAN/ACS/080	Instrument Cooling Fan - All Models	£39.38	
FAN/ACS/120	Cooling Fan -Compact Models 120mm <i>All C40 and H60 Models and C60 Models after S/No: 3786</i>	£43.94	
FAN/MID/C56	Cooling Fan - Midas 52/55/56 model	£84.26	
FAN/QCS/000	Cooling Fan - QCS and RSC Models <i>12in (305mm) Diameter</i>	£108.37	
FAN/TAN/C60	Cooling Fan - C60 model <i>Tangential Fan - Models with Serial No up to 3786</i>	£98.01	
FLO/WAT/SWT	Float Switch For Vac Pump Header Tank	£33.73	
FLS/ACS/316	Float Switch Ball Type S/S <i>RSC Water Level</i>	£334.12	
GAS/ACS/500	Door Seal Circular 100-200L Models	£119.25	
GAS/ACS/Q63	Door Seal Circular 320-400L Models	£129.02	
GAS/BOX/100	Door Seal Boxer 100L Blue	£130.97	
GAS/BOX/200	Door Seal Boxer 200L Blue	£152.82	















Part Number	Description	Price	
GAS/MID/300	Door Gasket - Midas V Series <i>Midas 32/36 & 55/56 models</i>	£162.89	
GAS/QCA/350	Door Seal Circular Compact Models	£82.56	
GAS/RSC/230	Door Seal Rectangular 230-350L Models	£161.33	
GAS/RSC/450	Door Seal Rectangular 450-700L Models	£183.22	
GAS/STG/001	Gasket -Steam Generator 680mm (Large) <i>Viton (fkm) 680mm ID joined rings For 'Barrel style' vessels - Pre 2013</i>	£92.59	
GAS/STG/002	Gasket-Steam Generator 377mm (Small) <i>Viton (fkm) 377mm ID Moulded rings</i>	£59.48	
GAT/Q63/SAF	Q63 Safety Gate	£81.66	
GAT/QCS/SAF	Safety Gate QCS	£59.60	
HBA/GLA/001	Hinge Bush Bearing RSC & Q63 <i>N.B. For Q63 models the bearing must be trimmed before use/despatch</i>	£15.96	
HBA/RSC/003	Hand Bolt Complete RSC Bottom (THERMAL)	£297.29	
HBA/WHE/REP	Hand-Wheel <i>Hand-Wheel only - Undrilled</i>	£84.65	
HBR/RSC/000	Threaded Insert Brass RSC	£80.71	
HCO/ACS/150	Thermal Fuse 140C	£6.01	
HCO/STG/180	Thermal Fuse 182C	£14.92	

Part Number	Description	Price	
HEX/CON/015	Heat exchanger 1/2in BSP <i>Working pressure max to 14.3 bar</i>	£447.11	
HEX/CON/022	Heat exchanger 3/4in BSP <i>Working pressure max to 14.3 bar</i>	£635.62	
HEX/CON/028	Heat exchanger 1in BSP <i>Working pressure max to 14.3 bar</i>	£934.43	
HND/QCS/YEL	Door Handle - QCS & Compact Series	£43.69	
HOS/STM/015	Flexible Hose 1/2in Braided S/S 1500mm	£122.89	
HOS/STM/022	Flexible Hose 3/4in Braided S/S 1500mm	£170.45	
HOS/STM/028	Flexible Hose 1in Braided S/S 1500mm	£178.40	
HTR/ACS/HOR	Heater - EH Models 3.5KW 230V	£134.15	
HTR/ACS/VER	Heater - EV Models 3.5KW 230V	£134.15	
HTR/ISG/105	Heater - Steam Generator 10.5KW 400V <i>Heater 10.5kW@400V IPR14658/1 2x12kW steps. 6 elements Made by INCOLOY800, flange 145 mm x 12 mm with a threaded hole 1/8" made by AISI304</i>	£403.70	
HTR/ISG/240	Heater - Steam Gen 24 KW 400V <i>Heater 24kW@415V IPR14658/1 2x12kW steps. 6 elements Made by INCOLOY800, flange 145 mm x 12 mm with a threaded hole 1/8" made by AISI304</i>	£556.66	
HTR/ISG/480	Heater - Steam Generator 48 KW 400V	£787.54	
HTR/MAT/001	Heater Mat 350W 230V - Rectangular <i>UL Approved</i>	£75.53	
HTR/MAT/002	Heater Mat 350W 230V - Square <i>UL Approved</i>	£75.53	

Part Number	Description	Price	
HTR/MID/C40	Heater Benchtop & 320-400L Models 2.9kW 230V 13A	£138.18	
HTR/MID/C60	Heater Top Loading Compact Models 2.9kW 230V 13A	£138.18	
HYD/PLC/INT	PLC Interface Unit	£206.70	
INL/ACS/FIL	Inline Electrical Filter	£34.04	
KEY/SLK/001	Setting Lock Key Switch	£90.65	
KEY/SLK/003	Setting Lock Key ONLY	£33.74	
KEY/STP/455	Key - Emergency Stop	£19.26	
KEY/TLK/002	Thermal Lock Key Switch - TACTROL2	£47.07	
KEY/TLK/004	Thermal Lock Key ONLY - TACTROL2	£33.74	
KEY/TLK/005	Thermal Lock Switch Assembly <i>Includes Key-switch, Key, Mounting Bezel and contact Blocks</i>	£95.08	
LEV/ACS/OMR	Float-less Level Switch	£104.39	
LOO/QCS/DOR	Door Wiring Loom QCA/QCS	£224.26	
LWP/ACS/00C	Low Water Probe Type C all Models <i>All Models</i>	£155.33	
LWP/ACS/EXC	Low Water Probe Extended <i>For Steam Generators with extended tab.</i>	£155.33	








Part Number	Description	Price	
MCB/ACS/006	Circuit Breaker 6Amp 1 Pole Type B <i>MCB 10kA, B Curve 1 Pole, 6 Amp VDE</i>	£26.29	
MCB/ACS/016	Circuit Breaker 16A 3 Pole Type B <i>MCB 10kA, B Curve 3Pole 16A, 415VAC VDE</i>	£80.01	
MCB/ACS/032	Circuit Breaker 32A 3 Pole Type B <i>MCB 10kA, B Curve 3Pole 32A, 415VAC VDE</i>	£176.31	
MCB/ACS/040	Circuit Breaker 40A 3 Pole Type B	£176.31	
MCB/ACS/04C	Circuit Breaker 4Amp 3 Pole Type C Vac	£70.70	
MCB/ACS/306	Circuit Breaker 6Amp 3 Pole Type B <i>MCB 10kA, B Curve 3Pole 6A, 415VAC VDE</i>	£74.16	
MCB/ACS/63A	Circuit Breaker 63Amp 3 Pole Type B <i>MCB 10kA, B Curve 3Pole 63A, 415VAC</i>	£232.66	
MCB/ACS/C02	Circuit Breaker 2Amp 3 Pole Type C	£107.86	
MCB/ACS/C06	Circuit Breaker 6Amp 3 Pole Type C Vac	£79.49	
MCB/ACS/C08	Circuit Breaker 8Amp 3 Pole Type C Vac <i>Suitable for Vac Pump application</i>	£78.93	
MCB/ACS/C20	Circuit Breaker 20Amp 3 Pole Type B	£72.46	
MCB/VAC/001	Circuit Breaker 6Amp 1 Pole Type C <i>Suitable for Vacuum Pump application</i>	£40.09	
MCB/VAC/003	Circuit Breaker 10Amp 3 Pole Type C Vac <i>Suitable for Vac Pump application</i>	£78.46	
MPS/MID/001	Mains Power Switch/Circuit Breaker Green <i>Compact Autoclaves</i>	£55.07	

Part Number	Description	Price	
MSA/CAN/HOR	Microswitch Assembly Horizontal Canopy	£19.31	
MSA/CAN/VER	Microswitch Assembly Horizontal Vertical	£19.31	
MSA/DCL/QCS	Microswitch Cross Roller QCS Door Closed	£53.51	
MSA/LID/QCS	Microswitch Small <i>Top and Front Loaders Door/ Canopy/Bench-Top and Compact Door</i>	£13.04	
MSA/MID/SOL	Microswitch Solenoid Compact Models	£13.71	
MSA/QCK/SBT	Microswitch Waterproof Safety Bolt <i>Midas Quickseal Safety Bolt Microswitch</i>	£64.65	
MSA/RSC/ROL	Microswitch Roller Rectangular Models	£29.75	
MSA/RSV/COI	Microswitch Coil Spring Actuator RSV	£165.23	
MSA/RSV/MIC	Microswitch IP64 Power Door Models	£48.78	
MSA/RSV/PER	Microswitch IP30 Perpendicular_RSV	£32.84	
MSA/RSV/PER	Microswitch Perpendicular RSV Models	£39.23	
MSA/RSV/ROT	Microswitch Side Rotary Lever RSV Models	£60.85	
MSA/SBT/001	Microswitch Special Safety Bolt	£42.37	
MSA/SBT/002	Microswitch Special Safety Bolt <i>Safety Bolt - Special - Microswitch - Safety Bolt - Special</i>	£41.15	















Part Number	Description	Price	
MSA/SBT/ACS	Microswitch Large <i>Bench-Top, Compact & Rectangular Model Doors</i>	£30.32	
PBS/ACS/001	Switch Latching Push Button <i>Includes contact block</i>	£194.98	
PNA/PWS/100	Purified Water Level Sensor	£926.47	
PRB/ACS/003	Thermocouple Probe Duplex K Type	£99.22	
PRB/ACS/004	Thermocouple Probe Wandering K Type <i>Please note: Wandering probes are sold packs of 2</i>	£226.52	
PRB/ACS/005	Thermocouple Probe Load Simulator K Type	£50.87	
PRB/ACS/007	Thermocouple Probe Duplex 250mm K Type	£110.98	
PRB/ACS/FSS	Fitting Bulkhead Duplex ProbeT/C 1/8 S/S <i>For thermocouples</i>	£10.24	
PRB/PLU/SOC	Plug & Socket K Type Wandering Probe	£10.25	
PRB/PRT/003	Probe Duplex PT100 6 x70mm <i>6mm replaces 3mm probe from Serial No 4800. For first change on models before serial No. 4800 a new fitting (PRB/PRT/FIT) must be supplied</i>	£87.61	
PRB/PRT/004	PRT Probe Wandering	£113.26	
PRB/PRT/005	PRT Probe Load Simulator	£63.02	
PRB/PRT/006	PRT Probe Wandering Armoured	£143.78	
PRB/PRT/007	PRT Probe Simplex Length 250mm	£121.61	

Part Number	Description	Price	
PRB/SST/PRB	Grade 316 Stainless Steel Cable Tie Pk3 <i>For attaching heater protection probes to heaters. Size: 127 x 4.6 mm. Requires specialist fitting tool.</i>	£125.26	
PSG/ACS/100	Pressure Gauge 100mm (-1 to 3 bar)	£156.16	
PSG/ACS/B10	Pressure Gauge 100mm Boxer Autoclaves	£153.02	
PSG/GAU/063	Pressure Gauge 63mm (-1 to 3 bar) <i>Bench-Top and Compact Models</i>	£68.02	
PSG/STG/063	Pressure Gauge 63mm (0 - 5 bar)	£78.66	
PSS/ISG/LOW	Pressure Switch Controller LP RANCO <i>Pressure Controller for Steam Generators</i>	£62.74	
PSS/SWI/001	Pressure Switch Range 1.5 to 15 psi	£113.36	
PSS/SWI/FIT	Fitting Compression Str.1/8BSP Fim x1/8 <i>Pressure Switch Adaptor - 1/8 MPT to Capillary Pipe (Compression joint).</i>	£7.77	
PST/ACS/001	Pressure Transducer 0 to 4 Bar <i>Autoclave Pressure</i>	£174.82	
PST/ISG/001	Pressure Transducer 0 to 6 Bar <i>Steam Generators</i>	£174.82	
REG/AIR/012	Regulator Filter 1/2" manual drain <i>For Air Ballast</i>	£166.81	
RLY/008/12V	Relay 8 Pin 12VDC <i>100-320L Door Opening</i>	£17.79	
RLY/240/008	Relay 230V - 8 pin	£38.90	
RLY/240/011	Relay 240V - 11 Pin	£38.90	


Part Number	Description	Price	
RLY/MID/HTR	Relay Heater Control - Compact Series	£38.90	
SBG/ACS/015	Silicon Bung 15mm	£3.05	
SBG/ACS/022	Silicon Bung 22mm	£5.34	
SHM/COM/SET	Shim Set Bench-Top & Compact Models <i>Contents: 5 of each shim, Hinge and Clamp 1mm and 0.5mm</i>	£71.58	
SHM/Q63/SET	Shim Set 320-400L Models <i>Contents: 5 of each shim, Hinge and Clamp 1mm and 0.5mm</i>	£71.58	
SHM/QCS/SET	Shim Set 100-200L Models <i>Contents: 5 of each shim, Hinge and Clamp 1mm and 0.5mm</i>	£71.58	
SHM/RSC/001	Shim Hinge Block 230-700L Models 1mm <i>Stainless Steel</i>	£3.58	
SHM/RSC/005	Shim Hinge Block 230-700L Models 0.5mm <i>Stainless Steel</i>	£3.48	
SOL/MID/001	12V Door Locking Solenoid Kit. 40-85L <i>Replaces 240V Solenoid on Compact Models before Serial Number 4438. Includes 12V Solenoid, Actuators, Bracket, Microswitch and 12V power supply</i>	£224.20	
SOL/QCS/001	12V Door Locking Solenoid Kit.100-400L <i>Replaces Door Solenoid on QCS and Q63 Models before Serial Number 4438. Includes 12V Solenoid, Actuator and Bracket, If replacing a 240V solenoid a power supply (TAC/POW/002) is required and must be ordered separately.</i>	£149.13	
SOL/QCS/002	Door Lock Solenoid 12V 40-400L <i>For Top and Front Loading, Bench-Top and Compact Models - Solenoid Only. Serial Nos from 4438.</i>	£125.67	
SOL/RSC/001	Solenoid Door Locking RSC <i>Obsolete. For models before March 2020 S/No: 4787 Order Part No: SOL/RSC/002 with 24V Power Supply TAC/POW/003 to replace rectifier in control box</i>	£323.45	
SOL/RSC/002	Solenoid Door Locking 24V RSC <i>Standard on All models after March 2020 S/No 4787.</i>	£314.03	
SPG/LID/C60	Spring Door Lift C60 Model	£77.13	

Part Number	Description	Price	
SPG/MID/SOL	Solenoid Spring Return Compact Models	£7.35	
SPG/QCS/LRG	Spring Large DoorLift Top Loading Models	£79.55	
STP/ACS/001	Emergency Stop Button <i>Without Contact Blocks</i>	£83.06	
STP/ACS/002	Emergency Stop Contact Block	£42.16	
STP/QCS/VER	Door Stop Rubber Black Vertical	£7.17	
STR/ACS/015	Strainer Brass 15mm	£53.88	
STR/ACS/022	Strainer Brass 22mm	£86.05	
STR/ACS/B22	Mesh Basket 3/4" Drain Filters <i>For threaded fittings only</i>	£27.36	
STR/ACS/S15	Strainer Inline S/S 15mm	£83.26	
STR/ACS/S22	Strainer Inline S/S 22mm	£103.84	
STR/ACS/WHL	Conical Strainer Large <i>Fits 28mm Tank Connectors</i>	£31.10	
STR/ACS/WHM	Conical Strainer Medium <i>Fits 22mm Tank Connectors</i>	£27.36	
STR/ACS/WHS	Conical Strainer Small <i>Fits 15mm Tank Connectors</i>	£24.19	
T931625PRT	Probe Duplex PT100 6mm x 250mm <i>6mm replaces 3mm probe from Serial No 4800. For first change on models before serial No. 4800 a new fitting (PRB/PRT/FIT) must be supplied</i>	£87.61	













Part Number	Description	Price	
TAC/OPT/001	Tactrol Output Board	£218.21	
TAC/OPT/FUS	Fuse 3.15A 5x20mm Tactrol O/P Board	£6.59	
TAC/PCB/001	Tactrol Main Board	£1,533.06	
TAC/PCB/021	Main Board Tactrol 21(IOT) <i>For Tactrol 21 models</i>	£1,643.45	
TAC/PCB/21M	Mini Board Tactrol21 IOT-MINI	£477.00	
TAC/PCB/MIN	Tactrol Mini Processor Main Board	£566.70	
TAC/PNL/001	Tactrol 2 Panel Printer Blank	£149.13	
TAC/POW/001	Power Supply Switch Mode	£243.90	
TAC/POW/002	Power Supply 240V-12V Solenoid Door	£50.12	
TAC/POW/003	Power supply 240V -24V 60W RSC Door <i>To drive replacement 24V Door Lock Solenoid for Rectangular Models. Required with Solenoid for models before Serial No. 4787 - March 2020.</i>	£54.82	
TAC/PRN/000	Tactrol Thermal Printer <i>Includes Thermal Printer, Serial Interface Adaptor and Mounting Plate</i>	£850.66	
TAC/PRN/003	Printer Type 3 24v Boxer Models	£829.92	
TAC/PRN/004	PrinterTactrol Touch Thermal <i>For direct connection to Tactrol Touch PLC System</i>	£850.66	
TAC/PRN/TUP	Thermal Printer Upgrade Kit <i>Includes Thermal Printer (TAC/PRN/002) plus associated Communication Adapter and Mounting Panel</i>	£477.21	

Part Number	Description	Price	
TAC/RES/BTN	Tactrol Reset Button	£10.56	
TAC/RIB/DIS	Ribbon Cable Display to Main Board Long <i>Models with Tactrol Board Not Mounted in Control Box</i>	£68.37	
TAC/RIB/DISS	Ribbon Cable Display to Main Board Short <i>Models with Tactrol Board Mounted in Control Box</i>	£25.94	
TAC/RIB/PRIL	Ribbon Cable Impact Printer Long <i>Models without a control box</i>	£49.14	
TAC/RIB/PRIS	Ribbon Cable Impact Printer Short <i>Models with a control box</i>	£27.31	
TAC/RIB/PRN	Ribbon Cable Printer to Main Board <i>Power Door and Double Door Models. Serial Number required - length varies with model</i>	£38.22	
TAC/RIB/PRNL	Ribbon Cable Thermal Printer Long <i>Models without a control box - Thermal Printer</i>	£49.14	
TAC/RIB/PRNS	Ribbon Cable Thermal Printer Short <i>Models with a control box</i>	£27.31	
TAC/USB/001	Tactrol USB Interface Card	£242.94	
TAC/USB/COV	Tactrol USB Black Dust Cap	£7.08	
TC2/CON/001	Tactrol 2 Main Control Panel	£486.02	
TC2/PRG/005	Tactrol 2 Programmer Panel	£590.29	
TC3/CON/001	Tactrol 3 Main Control Panel	£894.79	
THW/ACS/002	Load Simulator Thermowell Compact /QCS	£89.84	

Part Number	Description	Price	
THW/MID/001	Load Simulator Thermowell Compact Models	£84.69	
TIM/ACS/DIG	Digital Timer	£1,076.59	
TMP/ACS/001	Temperature Controller <i>Jumo Itron8</i>	£825.52	
TRA/480/230	Transformer 480V 230V / 3.5A	£1,456.49	
TST/STG/210	Thermostat 40 to 210C <i>Drain Condensers and Steam Generators</i>	£37.58	
VAC/DIA/DBL	Vacuum Pump Twin Head Diaphragm <i>40-60L Models</i>	£1,211.65	
VAC/DIA/SGL	Vacuum Pump Single Head Diaphragm	£1,048.17	
VAC/LRP/LRG	Vacuum Pump Liquid Ring 3 phase 450-700L <i>Q63 or RSC Models</i>	£2,051.25	
VAC/LRP/MED	Vacuum Pump Liquid Ring 3 phase 230-400L	£1,824.50	
VAC/LRP/SM1	Vacuum Pump Liquid Ring 1 phase 100-200L <i>100-200L Models</i>	£1,588.70	
VAC/LRP/SM3	Vacuum Pump Liquid Ring 3 phase 100-200L <i>100-200L Models</i>	£1,508.56	
VCE/ACS/001	Condensate Trap BPT13S 1/2in 15mm	£196.65	
VCE/ACS/002	Element Steam Condensate Trap (E-Fill)	£171.51	
VFL/ACS/015	Flapper Valve Compression 15mm	£17.98	

Part Number	Description	Price	
VFL/ACS/KIT	Hydroflow replacement valve fitting kit <i>For serial numbers before 5188 11-11-22. Adaptor kit for correct fitting to vacuum pump header tanks. Consists of 1x 15mm FIC fitting (ADT/ACS/15F), 1x Roll of PTFE tape (ZZZ/TAP/001) and an instruction sheet.</i>	£51.57	
VFL/ACS/TOR	Valve Hydroflo/Torbeck Side Inlet (NEW) <i>From S/No: 5188 11-11-22. Replacements for serial numbers before this serial will need an adaptor kit (VFL/ACS/KIT) for correct fitting to vacuum pump header tanks.</i>	£31.43	
VGA/ACS/015	Valve Gate 15mm	£38.53	
VGA/ACS/022	Valve Gate 22mm	£63.64	
VHA/090/015	Valve 15mm 1/4 Turn Ball 15mm	£37.04	
VHA/090/022	Valve 15mm 1/4 Turn Ball 22mm	£64.43	
VNR/ACS/015	Valve Non-Return 15mm (1/2in)	£24.84	
VNR/ACS/022	Valve Non-Return 22mm (3/4in)	£30.04	
VNR/ACS/028	Valve Non-Return 28mm (1in)	£34.25	
VPC/ACS/001	Air Purge Valve	£265.08	
VPE/ACS/002	Element Air Purge Valve (J-fill)	£219.60	
VSC/ACS/015	Safety Valve 15mm (1/2in) 2.5Bar	£231.83	
VSC/ACS/022	Safety Valve 22mm (3/4in) 2.5Bar	£302.34	

Part Number	Description	Price	
VSC/ASM/015	Safety Valve 15mm, 2.4Bar US Models (UL) <i>Approval acc. to ASME Boiler and Pressure Vessel Code, Section VIII.Div 1 (ASME)</i>	£265.47	
VSC/SSM/015	Safety Valve 1/2" x1" Stainless ASME <i>Approval acc. to ASME Boiler and Pressure Vessel Code, Section VIII.Div 1 (ASME)</i>	£582.30	
VSC/STG/022	Safety Valve 4.3bar 22mm Steam Generator	£296.62	
VVC/ACS/015	Valve Solenoid 15mm Vent	£134.63	
VVC/ACS/022	Valve Solenoid 22mm	£141.34	
VVC/PNU/015	Valve Pneumatic 15mm (1/2in)	£182.29	
VVC/PNU/022	Valve Pneumatic 22mm (3/4in)	£223.72	
VVC/PNU/028	Valve Pneumatic 28mm (1in)	£270.97	
VVC/RSC/028	Valve Solenoid 28mm <i>For serial numbers lower than 4628 please order the conversion kit VVC/RSC/ADP to ensure the correct valve length for direct replacement.</i>	£223.44	
VVC/RSC/ADP	28mm Solenoid Valve Converter Kit <i>Adapter Kit to fit 28mm Solenoid Valve to Autoclaves before serial number 4628.</i>	£30.93	
VVF/ACS/015	Fitting Vent Valve 15mm x 1/2 MIC	£4.36	
VVF/ACS/022	Fitting Vent Valve 22mm x 3/4 MIC	£4.94	
VVF/ACS/028	Fitting Vent Valve 28mm x 1in MIC	£8.05	
WAS/FIB/012	Washer Fibre - Heaters	£1.81	

Part Number	Description	Price	
WAS/PAC/015	Fibre Washers 15mm Pack of 10	£10.28	
WAS/PAC/022	Fibre Washers 22mm Pack of 10	£12.19	
WAS/PAC/028	Fibre Washers 28mm Pack of 10	£14.63	
WAS/SET/001	Fibre Washer Set <i>Complete set of washers, with PTFE tape & jointing compound</i>	£58.46	
WAT/ISG/PMP	Water Pump - Steam Generators	£377.51	
ZZZ/HYL/001	Hylomar Tube	£32.52	
ZZZ/TAP/001	Tape PTFE 12Mx12mm roll	£6.94	
ZZZ/TCG/001	Grease Copper Tube	£30.00	
ZZZ/TSG/001	Grease Silicone Tube 100g	£32.52	
ZZZ/TSS/001	Silicon Sealant Translucent / Clear	£19.55	
ZZZ/TSS/002	Silicon Sealant Black	£19.55	
ZZZ/TSS/003	Silicon Sealant White	£19.55	

Part Number	Description	Price	
Consumable Items			
ZZZ/BDT/000	Bowie Dick Test Pack BD-1XE Pack of 1	£25.00	
ZZZ/DEO/001	Deodorising Capsules -Classic/Peach <i>Please contact service/sales ordering this two.</i>	£95.00	
ZZZ/DEO/002	Deodorising Capsules -Apple	£95.00	
ZZZ/DEO/003	Deodorising Capsules -Citrus	£95.00	
ZZZ/LOG/000	Priorclave Autoclave Log Book	£25.00	
ZZZ/MSS/06E	Biological Indicator MSS-06E Pack of 100	£300.00	
ZZZ/RIB/002	Tactrol Printer Ribbon Type 2	£23.45	
ZZZ/ROL/005	Printer Roll (Pk. 5)	£23.45	
ZZZ/ROL/007	Thermal Printer Roll (Pk 5)	£28.19	
ZZZ/SYP/LGE	Siphon Pump Large <i>Long Siphon Pump for Top Loading Autoclaves</i>	£50.00	
ZZZ/SYP/SML	Siphon Pump Small <i>Short Siphon Pump for Front Loading Autoclaves</i>	£30.00	
ZZZ/TST/S-E	Indicator TST Class 6 Autoclave S-E <i>Pack of 250 Test Strips</i>	£55.00	

Please Note

Prices shown are List, EX-Works, UK and do not include packing or shipping costs, which will be added as appropriate or any duties payable outside of the UK.