




























# Illustrated Spares Price List

















18/01/2021




Part Number	Description	Price	
<b>AFL/ACS/EXH</b>	Containment Filter Element - External <i>For Filter Housing</i>	£431.63	
<b>AFL/ACS/EXI</b>	Containment Filter Element - Internal	£550.69	
<b>AFL/ACS/SML</b>	Filter - Air Intake - SMALL	£129.52	
<b>CAS/MID/000</b>	Castor Wheel Set - Compact C60 model <i>Set of 4</i>	£27.61	
<b>CON/230/24A</b>	Contactor - 50amp	£92.99	
<b>CON/230/74A</b>	Contactor-74 Amp 230V	£250.98	
<b>CON/240/09A</b>	Contactor - 9 Amp 240V	£45.68	
<b>CON/240/23A</b>	Contactor - 23 Amp 240V	£78.51	
<b>CON/240/45A</b>	Contactor - 100 Amp 240V	£173.95	
<b>CON/240/811</b>	Overload Relay (8A-11A Range)	£58.94	
<b>CON/24V/09A</b>	Contactor - 9 Amp 24V	£57.25	
<b>CON/OLD/004</b>	Overload Relay (2.4 - 4.0 Amp range)	£40.40	
<b>CON/OLD/69A</b>	Overload Relay (6-9 Amp range)	£71.80	











Part Number	Description	Price	
CON/OLD/OL6	Overload Relay (4-6 Amp range)	£33.54	
CTR/ACS/001	Cycle Counter Panel Mount 12V DC	£29.94	
DTM/240/002	Delay Timer - 2 Contact	£303.78	
FAN/ACS/080	Instrument Cooling Fan - All Models	£36.07	
FAN/ACS/120	Cooling Fan -Compact Models 120mm <i>All C40 and H60 Models and C60 Models after S/No: 3786</i>	£40.24	
FAN/MID/C56	Cooling Fan - Midas 52/55/56 model	£77.17	
FAN/QCS/000	Cooling Fan - QCS and RSC Models <i>12in (305mm) Diameter</i>	£97.83	
FAN/TAN/C60	Cooling Fan - C60 model <i>Tangential Fan - Models with Serial No up to 3786</i>	£89.77	
FLO/WAT/SWT	Float Switch - Vac Pump Header Tank	£30.89	
FLS/ACS/316	Float Switch Ball Type S/S <i>RSC Water Level</i>	£306.03	
GAS/ACS/500	Door Seal Circular 100-200L Models	£109.22	
GAS/ACS/Q63	Door Seal Circular 320-400L Models	£118.17	
GAS/MID/300	Door Gasket - Midas V Series <i>Midas 32/36 &amp; 55/56 models</i>	£149.19	
GAS/QCA/350	Door Gasket QCA Compact Series	£75.62	















Part Number	Description	Price	
GAS/RSC/230	Door Seal Rectangular 230-350L Models	£147.77	
GAS/RSC/450	Door Seal Rectangular 450-700L Models	£167.82	
GAS/STG/001	Gasket -Steam Generator 680mm (Large) <i>Viton (fkm) 680mm ID joined rings For 'Barrel style' vessels - Pre 2013</i>	£84.81	
GAS/STG/002	Gasket-Steam Generator 377mm (Small) <i>Viton (fkm) 377mm ID Moulded rings</i>	£54.48	
GAT/Q63/SAF	Q63 Safety Gate	£74.80	
GAT/QCS/SAF	Safety Gate QCS	£54.59	
HBA/GLA/001	Hinge Bush Bearing RSC & Q63 <i>N.B. For Q63 models the bearing must be trimmed before use/despatch</i>	£13.82	
HBA/RSC/001	Hand Bolt Complete RSC Top (SAFETY)	£280.61	
HBA/RSC/003	Hand Bolt Complete RSC Bottom (THERMAL)	£272.29	
HBA/WHE/REP	Hand-Wheel <i>Hand-Wheel only - Undrilled</i>	£77.53	
HBR/ACS/002	Threaded Insert Brass RSC	£73.92	
HCO/ACS/150	Thermal Fuse 150°C	£5.50	
HCO/STG/180	Thermal Fuse 182C	£13.67	
HEX/CON/015	Heat exchanger 1/2in BSP	£409.51	

Part Number	Description	Price	
HEX/CON/022	Heat exchanger 3/4in BSP	£582.17	
HEX/CON/028	Heat exchanger 1in BSP	£855.86	
HND/QCS/YEL	Door Handle - QCS & Compact Series	£40.02	
HOS/STM/015	Flexible Hose 1/2in Braided S/S 1500mm	£112.55	
HOS/STM/022	Flexible Hose 3/4in Braided S/S 1500mm	£150.49	
HOS/STM/028	Flexible Hose 1in Braided S/S 1500mm	£163.40	
HTR/ACS/HOR	Heater - EH Models 3.5KW 230V	£122.87	
HTR/ACS/VER	Heater - EV Models 3.5KW 230V	£122.87	
HTR/ISG/105	Heater - Steam Generator 10.5KW 400V <i>Heater 10.5kW@400V IPR14658/1 2x12kW steps. 6 elements Made by INCOLOY800, flange 145 mm x 12 mm with a threaded hole 1/8" made by AISI304</i>	£310.34	
HTR/ISG/240	Heater - Steam Gen 24 KW 400V <i>Heater 24kW@415V IPR14658/1 2x12kW steps. 6 elements Made by INCOLOY800, flange 145 mm x 12 mm with a threaded hole 1/8" made by AISI304</i>	£509.85	
HTR/ISG/480	Heater - Steam Generator 48 KW 400V	£721.32	
HTR/MAT/001	Heater Mat 350W 230V - Rectangular <i>UL Approved</i>	£69.17	
HTR/MAT/002	Heater Mat 350W 230V - Square <i>UL Approved</i>	£69.17	
HTR/MID/C40	Heater Benchtop & 320-400L Models 2.9kW 230V 13A	£126.56	

Part Number	Description	Price	
HTR/MID/C60	Heater Top Loading Compact Models 2.9kW 230V 13A	£126.56	
INL/ACS/FIL	Inline Electrical Filter	£31.17	
KEY/SLK/001	Setting Lock Key Switch	£83.03	
KEY/SLK/003	Setting Lock Key ONLY	£30.90	
KEY/STP/455	Key - Emergency Stop	£17.64	
KEY/TLK/002	Thermal Lock Key Switch - TACTROL2	£43.12	
KEY/TLK/004	Thermal Lock Key ONLY - TACTROL2	£30.90	
KEY/TLK/005	Thermal Lock Switch Assembly <i>Includes Key-switch, Key, Mounting Bezel and contact Blocks</i>	£87.09	
LEV/ACS/OMR	Float-less Level Switch	£95.61	
LOO/QCS/DOR	Door Wiring Loom QCA/QCS	£205.41	
LWP/ACS/00C	Low Water Probe Type C all Models <i>All Models</i>	£142.27	
LWP/ACS/EXC	Low Water Probe Extended <i>For Steam Generators with extended tab.</i>	£142.27	
LWP/DPB/000	Low Water Ultrasonic Sensor <i>Dual low water probe for boilers. Pure water supplies</i>	£848.57	
MCB/ACS/006	Circuit Breaker 6Amp 1 Pole Type B	£24.08	















Part Number	Description	Price	
MCB/ACS/016	Circuit Breaker 16A 3 Pole Type B	£73.28	
MCB/ACS/032	Circuit Breaker 32A 3 Pole Type B	£161.49	
MCB/ACS/040	Circuit Breaker 40A 3 Pole Type B	£161.49	
MCB/ACS/045	Circuit Breaker Type B50 3-Pole	£203.88	
MCB/ACS/04C	Circuit Breaker 4Amp 3 Pole Type C Vac	£64.76	
MCB/ACS/306	Circuit Breaker 6Amp 3 Pole Type B	£67.92	
MCB/ACS/63A	Circuit Breaker 63Amp 3 Pole Type B	£213.10	
MCB/ACS/C02	Circuit Breaker 2Amp 3 Pole Type C	£98.79	
MCB/ACS/C06	Circuit Breaker 6Amp 3 Pole Type C Vac	£72.81	
MCB/ACS/C08	Circuit Breaker 8Amp 3 Pole Type C Vac <i>Suitable for Vac Pump application</i>	£72.29	
MCB/ACS/C20	Circuit Breaker 20Amp 3 Pole Type B	£66.37	
MCB/VAC/001	Circuit Breaker 6Amp 1 Pole Type C <i>Suitable for Vacuum Pump application</i>	£36.72	
MCB/VAC/003	Circuit Breaker 10Amp 3 Pole Type C Vac <i>Suitable for Vac Pump application</i>	£71.86	
MPS/MID/001	Mains Power Switch/Circuit Breaker Green <i>Compact Autoclaves</i>	£50.44	




Part Number	Description	Price	
MSA/CAN/HOR	Microswitch Assembly Horizontal Canopy	£17.69	
MSA/CAN/VER	Microswitch Assembly Horizontal Vertical	£17.69	
MSA/DCL/QCS	Microswitch Cross Roller QCS Door Closed	£49.01	
MSA/LID/QCS	Microswitch Small <i>Top and Front Loaders Door/ Canopy/Bench-Top and Compact Door</i>	£11.94	
MSA/MID/SOL	Microswitch Solenoid Compact Models	£11.94	
MSA/QCK/SBT	Microswitch Waterproof Safety Bolt <i>Midas Quickseal Safety Bolt Microswitch</i>	£59.21	
MSA/RSC/ROL	Microswitch Roller Rectangular Models	£27.25	
MSA/RSV/COI	Microswitch Coil Spring Actuator RSV	£73.43	
MSA/RSV/MIC	Microswitch IP64 Power Door Models	£44.68	
MSA/RSV/PER	Microswitch Perpendicular RSV Models	£35.93	
MSA/RSV/ROT	Microswitch Side Rotary Lever RSV Models	£55.74	
MSA/SBT/001	Microswitch Special Safety Bolt	£38.81	
MSA/SBT/002	Microswitch Special Safety Bolt <i>Safety Bolt - Special - Microswitch - Safety Bolt - Special</i>	£37.69	
MSA/SBT/ACS	Microswitch Large <i>Bench-Top, Compact &amp; Rectangular Model Doors</i>	£27.77	

Part Number	Description	Price	
PBS/ACS/001	Switch Latching Push Button <i>Includes contact block</i>	£178.58	
PNA/PWS/100	Purified Water Level Sensor	£848.57	
PRB/ACS/003	Thermocouple Probe Duplex K Type	£90.87	
PRB/ACS/004	Thermocouple Probe Wandering K Type <i>Please note: Wandering probes are sold packs of 2</i>	£207.48	
PRB/ACS/005	Thermocouple Probe Load Simulator K Type	£46.59	
PRB/ACS/007	Thermocouple Probe Duplex 250mm K Type	£101.65	
PRB/ACS/FSS	Fitting Bulkhead Duplex Probe 1/8 S/S <i>For thermocouples</i>	£9.38	
PRB/PLU/SOC	Plug & Socket K Type Wandering Probe	£9.39	
PRB/PRT/003	Probe Duplex PT100 6 x70mm <i>6mm replaces 3mm probe from Serial No 4800. For first change on models before serial No. 4800 a new fitting (PRB/PRT/FIT) must be supplied</i>	£80.25	
PRB/PRT/004	PRT Probe Wandering <i>Please note: Wandering probes are sold packs of 2</i>	£207.48	
PRB/PRT/005	PRT Probe Load Simulator	£57.72	
PRB/PRT/006	PRT Probe Wandering Armoured <i>Please note: Wandering probes are sold packs of 2</i>	£263.37	
PRB/PRT/007	Probe Duplex PT100 6mm x 250mm <i>6mm replaces 3mm probe from Serial No 4800. For first change on models before serial No. 4800 a new fitting (PRB/PRT/FIT) must be supplied</i>	£80.25	
PRB/SST/PRB	Grade 316 Stainless Steel Cable Tie Pk3 <i>For attaching heater protection probes to heaters. Size: 127 x 4.6 mm. Requires specialist fitting tool.</i>	£114.73	

















Part Number	Description	Price	
PSG/ACS/100	Pressure Gauge 100mm (-1 to 3 bar)	£143.03	
PSG/ACS/163	Pressure Gauge 100mm Boxer Autoclaves	£140.16	
PSG/GAU/063	Pressure Gauge 63mm (-1 to 3 bar) <i>Bench-Top and Compact Models</i>	£52.15	
PSG/STG/063	Pressure Gauge 63mm (0 - 5 bar)	£72.05	
PSS/ISG/LOW	Pressure Switch Controller LP RANCO <i>Pressure Controller for Steam Generators</i>	£57.47	
PSS/SWI/001	Pressure Switch Range 1.5 to 15 psi	£103.83	
PSS/SWI/FIT	Fitting Compression Str.1/8BSP Fim x1/8 <i>Pressure Switch Adaptor - 1/8 MPT to Capillary Pipe (Compression joint).</i>	£7.12	
PST/ACS/001	Pressure Transducer 0 to 4 Bar <i>Autoclave Pressure</i>	£160.12	
PST/ISG/001	Pressure Transducer 0 to 6 Bar <i>Steam Generators</i>	£160.12	
REG/AIR/012	Regulator Filter 1/2" manual drain <i>For Air Ballast</i>	£152.79	
RLY/008/12V	Relay 8 Pin 12V DC Assembly <i>For QCS Door Opening</i>	£16.29	
RLY/240/008	Relay 230V - 8 pin	£35.63	
RLY/240/011	Relay 240V - 11 Pin	£35.63	
RLY/MID/HTR	Relay Heater Control - Compact Series	£35.63	

Part Number	Description	Price	
<b>SBG/ACS/015</b>	Silicon Bung 15mm	£2.80	
<b>SBG/ACS/022</b>	Silicon Bung 22mm	£4.89	
<b>SHM/COM/SET</b>	Shim Set Bench-Top & Compact Models <i>Contents: 5 of each shim, Hinge and Clamp 1mm and 0.5mm</i>	£65.56	
<b>SHM/Q63/SET</b>	Shim Set 320-400L Models <i>Contents: 5 of each shim, Hinge and Clamp 1mm and 0.5mm</i>	£65.56	
<b>SHM/QCS/SET</b>	Shim Set 100-200L Models <i>Contents: 5 of each shim, Hinge and Clamp 1mm and 0.5mm</i>	£65.56	
<b>SHM/RSC/001</b>	Shim Hinge Block 320-700L Models 1mm <i>Stainless Steel</i>	£3.28	
<b>SHM/RSC/005</b>	Shim Hinge Block 320-700L Models 0.5mm <i>Stainless Steel</i>	£3.18	
<b>SOL/MID/001</b>	12V Door Locking Solenoid Kit. 40-85L <i>Replaces 240V Solenoid on Compact Models before Serial Number 4438. Includes 24V Solenoid, Actuators, Bracket, Microswitch and 12V power supply</i>	£205.35	
<b>SOL/QCS/001</b>	12V Door Locking Solenoid Kit.100-400L <i>Replaces Door Solenoid on QCS and Q63 Models before Serial Number 4438. Includes 12V Solenoid, Actuator and Bracket, If replacing a 240V solenoid a power supply (TAC/POW/002) is required and must be ordered separately.</i>	£115.11	
<b>SOL/QCS/002</b>	Door Lock Solenoid 12V 40-400L <i>For Top and Front Loading, Bench-Top and Compact Models - Solenoid Only. Serial Nos from 4438.</i>	£136.59	
<b>SOL/RSC/001</b>	Solenoid Door Locking RSC <i>Obsolete. For models before March 2020 S/No: 4787 Order Part No: SOL/RSC/002 with 24V Power Supply TAC/POW/003 to replace rectifier in control box</i>	£296.25	
<b>SPG/LID/C60</b>	Spring Door Lift C60 Model	£70.64	
<b>SPG/MID/SOL</b>	Solenoid Spring Return Compact Models	£6.73	
<b>SPG/QCS/LRG</b>	Spring Large DoorLift Top Loading Models	£72.86	













Part Number	Description	Price	
STP/ACS/001	Emergency Stop Button <i>Without Contact Blocks</i>	£76.08	
STP/ACS/002	Emergency Stop Contact Block	£38.61	
STP/QCS/VER	Door Stop Rubber Black Vertical	£6.56	
STR/ACS/015	Strainer Brass 15mm	£49.35	
STR/ACS/022	Strainer Brass 22mm	£78.82	
STR/ACS/B22	Mesh Basket 3/4" Drain Filters <i>For threaded fittings only</i>	£25.06	
STR/ACS/S15	Strainer Inline S/S 15mm	£76.26	
STR/ACS/S22	Strainer Inline S/S 22mm	£95.11	
STR/ACS/WHL	Conical Strainer Large <i>Fits 28mm Tank Connectors</i>	£28.49	
STR/ACS/WHM	Conical Strainer Medium <i>Fits 22mm Tank Connectors</i>	£25.06	
STR/ACS/WHs	Conical Strainer Small <i>Fits 15mm Tank Connectors</i>	£22.16	
TAC/OPT/001	Tactrol Output Board	£199.86	
TAC/OPT/FUS	Fuse 3.15A 5x20mm Tactrol O/P Board	£6.04	
TAC/PCB/001	Tactrol Main Board	£1,404.16	

Part Number	Description	Price	
TAC/PCB/MIN	Tactrol Mini Processor Main Board	£519.05	
TAC/PNL/001	Tactrol 2 Panel Printer Blank	£136.59	
TAC/POW/001	Power Supply Switch Mode	£223.39	
TAC/POW/002	Power Supply 240V-12V Solenoid Door	£45.90	
TAC/POW/003	Power supply 240V -24V 60W RSC Door <i>To drive replacement 24V Door Lock Solenoid for Rectangular Models. Required with Solenoid for models before Serial No. 4787 - March 2020.</i>	£50.21	
TAC/PRN/000	Tactrol Thermal Printer <i>Includes Thermal Printer, Serial Interface Adaptor and Mounting Plate</i>	£779.14	
TAC/PRN/003	Printer Type 3 24v Boxer Models	£760.14	
TAC/PRN/004	Tactrol Touch Thermal Printer <i>For direct connection to Tactrol Touch PLC System</i>	£779.14	
TAC/PRN/TUP	Thermal Printer Upgrade Kit <i>Includes Thermal Printer (TAC/PRN/002) plus associated Communication Adapter and Mounting Panel</i>	£437.09	
TAC/RES/BTN	Tactrol Reset Button	£9.67	
TAC/RIB/DIS	Ribbon Cable Display to Main Board Long <i>Models with Tactrol Board Not Mounted in Control Box - Serial Number required - length varies with model</i>	£47.22	
TAC/RIB/DISS	Ribbon Cable Display to Main Board Short <i>Models with Tactrol Board Mounted in Control Box</i>	£23.76	
TAC/USB/001	Tactrol USB Interface Card	£222.51	
TAC/USB/COV	Tactrol USB Black Dust Cap	£6.49	












Part Number	Description	Price	
TC2/CON/001	Tactrol 2 Main Control Panel	£445.16	
TC2/PRG/005	Tactrol 2 Programmer Panel	£540.66	
TC3/CON/001	Tactrol 3 Main Control Panel	£819.55	
THW/ACS/002	Load Simulator Thermowell Compact /QCS	£82.28	
THW/MID/001	Load Simulator Thermowell Compact Models	£77.57	
TIM/ACS/DIG	Digital Timer	£986.07	
TMP/ACS/001	Temperature Controller <i>Jumo Itron8</i>	£756.11	
TST/STG/210	Thermostat 40 to 210C <i>Drain Condensers and Steam Generators</i>	£34.42	
VAC/DIA/DBL	Vacuum Pump Twin Head Diaphragm	£1,109.78	
VAC/DIA/SGL	Vacuum Pump Single Head Diaphragm	£960.04	
VAC/LRP/LRG	Vacuum Pump Liquid Ring 3 phase 450-700L <i>Q63 or RSC Models</i>	£1,878.78	
VAC/LRP/MED	Vacuum Pump Liquid Ring 3 phase 230-400L <i>100-700L Models</i>	£1,671.10	
VAC/LRP/SM1	Vacuum Pump Liquid Ring 1 phase 100-200L <i>100-200L Models</i>	£1,455.12	
VAC/LRP/SM3	Vacuum Pump Liquid Ring 3 phase 100-200L <i>100-200L Models</i>	£1,381.72	

Part Number	Description	Price	
VCE/ACS/001	Condensate Trap BPT13S 1/2in 15mm	£168.99	
VCE/ACS/002	Element Steam Condensate Trap (E-Fill)	£157.09	
VFL/ACS/015	Flapper Valve Compression 15mm	£16.47	
VFL/ACS/TOR	Valve Torbeck Type Auto Water Fill Tank	£28.79	
VGA/ACS/015	Valve Gate 15mm	£35.29	
VGA/ACS/022	Valve Gate 22mm	£58.29	
VHA/090/015	Valve 15mm 1/4 Turn Ball 15mm	£33.92	
VHA/090/022	Valve 15mm 1/4 Turn Ball 22mm	£59.01	
VNR/ACS/015	Valve Non-Return 15mm (1/2in)	£22.75	
VNR/ACS/022	Valve Non-Return 22mm (3/4in)	£27.51	
VNR/ACS/028	Valve Non-Return 28mm (1in)	£31.37	
VPC/ACS/001	Air Purge Valve	£216.67	
VPE/ACS/002	Element Air Purge Valve (J-fill)	£201.14	
VSC/ACS/015	Safety Valve 15mm (1/2in) 2.5Bar	£212.34	

Part Number	Description	Price	
VSC/ACS/022	Safety Valve 22mm (3/4in) 2.5Bar	£276.92	
VSC/STG/022	Safety Valve 4.3bar 22mm Steam Generator	£271.68	
VSI/ACS/015	Steam Inlet Valve 15mm	£157.52	
VVC/ACS/015	Valve Solenoid 15mm Vent	£123.31	
VVC/ACS/022	Valve Solenoid 22mm	£129.46	
VVC/PNU/015	Valve Pneumatic 15mm (1/2in)	£166.96	
VVC/PNU/022	Valve Pneumatic 22mm (3/4in)	£204.91	
VVC/PNU/028	Valve Pneumatic 28mm (1in)	£248.18	
VVC/RSC/028	Valve Solenoid 28mm <i>For serial numbers lower than 4628 please order the conversion kit VVC/RSC/ADP to ensure the correct valve length for direct replacement.</i>	£204.65	
VVC/RSC/ADP	28mm Solenoid Valve Converter Kit <i>Adapter Kit to fit 28mm Solenoid Valve to Autoclaves before serial number 4628.</i>	£28.33	
VVF/ACS/015	Fitting Vent Valve 15mm x 1/2 MIC	£3.99	
VVF/ACS/022	Fitting Vent Valve 22mm x 3/4 MIC	£4.52	
VVF/ACS/028	Fitting Vent Valve 28mm x 1in MIC	£7.37	
WAS/FIB/012	Washer Fibre - Heaters	£1.66	

Part Number	Description	Price	
WAS/PAC/015	Fibre Washers 15mm Pack of 10	£0.94	
WAS/PAC/022	Fibre Washers 22mm Pack of 10	£0.94	
WAS/PAC/028	Fibre Washers 28mm Pack of 10	£0.94	
WAS/SET/001	Fibre Washer Set <i>Complete set of washers, with PTFE tape &amp; jointing compound</i>	£53.54	
WAT/ISG/PMP	Water Pump - Steam Generators	£345.77	
ZZZ/HYL/001	Hylomar Tube	£29.79	
ZZZ/TAP/001	Tape PTFE 12Mx12mm roll	£6.36	
ZZZ/TCG/001	Grease Copper Tube	£25.75	
ZZZ/TSG/001	Grease Silicone Tube 100g	£29.79	
ZZZ/TSS/001	Silicon Sealant Translucent / Clear	£17.90	
ZZZ/TSS/002	Silicon Sealant Black	£17.90	
ZZZ/TSS/003	Silicon Sealant White	£17.90	



Part Number	Description	Price	
<b>Consumable Items</b>			
ZZZ/DEO/001	Deodorising Capsules -Classic/Peach <i>Please contact service/sales ordering this two.</i>	£87.55	
ZZZ/DEO/002	Deodorising Capsules -Apple	£87.55	
ZZZ/DEO/003	Deodorising Capsules -Citrus	£87.55	
ZZZ/LOG/000	Priorclave Autoclave Log Book	£20.60	
ZZZ/MSS/06E	Biological Indicator MSS-06E Pack of 100	£272.95	
ZZZ/RIB/002	Tactrol Printer Ribbon Type 2	£21.49	
ZZZ/ROL/005	Printer Roll (Pk. 5)	£21.49	
ZZZ/ROL/007	Thermal Printer Roll (Pk 5)	£25.82	
ZZZ/SYP/LGE	Siphon Pump Large <i>Long Siphon Pump for Top Loading Autoclaves</i>	£45.90	
ZZZ/SYP/SML	Siphon Pump Small <i>Short Siphon Pump for Front Loading Autoclaves</i>	£27.32	
ZZZ/TST/S-E	Indicator TST Class 6 Autoclave S-E <i>Pack of 100 Test Strips</i>	£51.50	

**Please Note**

Prices shown are List, EX-Works, UK and do not include packing or shipping costs, which will be added as appropriate or any duties payable outside of the UK.