

Priorclave Pass-Through autoclaves, also known as double ended sterilisers, enable sealed laboratory environments to be maintained, door interlocks preventing potential for laboratory standards to be compromised.

Every double-ended autoclave provides a sterile path in and out of sealed laboratories.

Each double-ended Priorclave is customer built at the company's UK manufacturing plant and is available in a range of chamber capacities and styles.

## KEY FEATURES

- Front loading Round and Rectangular chambers.
- Chamber capacities from 150 litres.
- Doors interlock to prevent opening both sides.
- TACTROL® microprocessor control system.
- Control panel both sides.
- Operation up to 138°C – 2.4bar.
- Forced air cooling helps reduce cycle times.
- Automatic timed freesteaming improves air removal and sterilising performance.
- Delayed start and media warming function.
- Epoxy coated panels and frame members are treated with an anti-bacterial agent which is effective against all bacteria and fungi including MRSA.
- Pressure vessel carries a 10 year warranty and is insurance approved.
- Service and support direct from manufacturer (or appointed distributor).

## FRONT LOADING PASS-THROUGH



## SPEAK WITH AN AUTOCLAVE EXPERT



+44 (0)20 8316 6620 • [priorclave.co.uk](http://priorclave.co.uk) • [sales@priorclave.co.uk](mailto:sales@priorclave.co.uk)

Our highly trained service technicians provide dedicated pre- and after-sales support including planned maintenance programmes to ensure you get the most from your Priorclave.