






EV100 OUTSIDE DIMENSIONS




EV100 WITH 1 LITER FLASKS




CALCULATE:

Approximately 16
Flask Bottles (1 L) size will
fit into 2 Baskets.


Wire Basket Dimensions
17.7" Dia x 8.3" D
450mm Dia x 210mm D
(2 stackable baskets provided)



VESSEL
DIMENSIONS
19.5" Dia x 19.5" D




EV100 WITH 2 LITER FLASKS




CALCULATE:

Approximately 4
Flask Bottles (2 L) size will
fit into 1 Basket.


Wire Basket Dimensions
17.7" Dia x 8.3" D
450mm Dia x 210mm D
(2 stackable baskets provided)



VESSEL
DIMENSIONS
19.5" Dia x 19.5" D




EV100 WITH 3 LITER FLASKS




CALCULATE:

Approximately 3
Flask Bottles (3 L) size will
fit into 1 Basket.


Wire Basket Dimensions
17.7" Dia x 8.3" D
450mm Dia x 210mm D
(2 stackable baskets provided)



VESSEL
DIMENSIONS
19.5" Dia x 19.5" D




EV100 WITH 4 LITER FLASKS




CALCULATE:

Approximately 3
Flask Bottles (4 L) size will
fit into 1 Basket.


Wire Basket Dimensions
17.7" Dia x 8.3" D
450mm Dia x 210mm D
(2 stackable baskets provided)



VESSEL
DIMENSIONS
19.5" Dia x 19.5" D




EV100 WITH 5 LITER FLASKS




CALCULATE:

Approximately 2
Flask Bottles (5 L) size will
fit into 1 Basket.


Wire Basket Dimensions
17.7" Dia x 8.3" D
450mm Dia x 210mm D
(2 stackable baskets provided)



VESSEL
DIMENSIONS
19.5" Dia x 19.5" D




EV100 WITH 6 LITER FLASKS




CALCULATE:

Approximately 1
Flask Bottles (6 L) size will
fit into 1 Basket.


Wire Basket Dimensions
17.7" Dia x 8.3" D
450mm Dia x 210mm D
(2 stackable baskets provided)



VESSEL
DIMENSIONS
19.5" Dia x 19.5" D




EV100 WITH 125 mL FLASKS




CALCULATE:

Approximately 66
Flask Bottles (125mL) size will
fit into 2 Baskets.


Wire Basket Dimensions
17.7" Dia x 8.3" D
450mm Dia x 210mm D
(2 stackable baskets provided)



VESSEL
DIMENSIONS
19.5" Dia x 19.5" D




EV100 WITH 250 mL FLASKS




CALCULATE:
Approximately 38
Flask Bottles (250mL) size will
fit into 2 Baskets.


Wire Basket Dimensions
17.7" Dia x 8.3" D
450mm Dia x 210mm D
(2 stackable baskets provided)



VESSEL
DIMENSIONS
19.5" Dia x 19.5" D



EV100 WITH 500 mL FLASKS



CALCULATE:


Approximately 26
Flask Bottles (500mL) size will
fit into 2 Baskets.

Wire Basket Dimensions


17.7" Dia x 8.3" D

450mm Dia x 210mm D

(2 stackable baskets provided)



VESSEL
DIMENSIONS
19.5" Dia x 19.5" D



EV100 WITH TEST TUBE RACKS

Typical Rack Size
Holds 60 Tubes



Tube size
18 - 20 mm

3.25"
82.6 mm


CALCULATE:

Approximately 8 Test
Tube Racks will fit
into 2 Baskets.


6.125"
155.6mm

9.625"
244.5mm

Wire Basket Dimensions
17.7" Dia x 8.3" D
450mm Dia x 210mm D
(2 stackable baskets provided)



VESSEL
DIMENSIONS
19.5" Dia x 19.5" D



EV100 WITH AUTOCLAVE
WASTE BAGS AND HOLDERS

Bag size measures 8.5" x 11.8"
215 x 300mm

Holder
measures
91 x 128
x 214mm

8.4"
214mm




CALCULATE:
Approximately 16
Waste Bags and
Holders will fit
into 2 Baskets.


5.0"
128mm

3.6"
91mm


Wire Basket Dimensions
17.7" Dia x 8.3" D
450mm Dia x 210mm D
(2 stackable baskets provided)



VESSEL
DIMENSIONS
19.5" Dia x 19.5" D




EV100 WITH AUTOCLAVE WASTE BAGS




CALCULATE:

Approximately 2 Large
Waste Bags will fit into
2 Baskets.



Wire Basket Dimensions
17.7" Dia x 8.3" D
450mm Dia x 210mm D
(2 stackable baskets provided)



VESSEL
DIMENSIONS
19.5" Dia x 19.5" D