



COMPACT 40 OUTSIDE DIMENSIONS




COMPACT 40 WITH 1 LITER FLASKS



CALCULATE:



Approximately 5
Flask Bottles (1 L) size will
fit on Bottom shelf.



BOTTOM SHELF DIMENSIONS


9.7" W x 16.5" D

COMPACT 40 WITH 2 LITER FLASKS




BOTTOM SHELF DIMENSIONS
9.7" W x 16.5" D

COMPACT 40 WITH 125mL FLASKS




CALCULATE:
Approximately 42
Flask Bottles (125mL) size will
fit on Top and Bottom shelves.




TOP SHELF DIMENSIONS
13.6"W x 16.5"D
WITH 6.9" CLEARANCE

BOTTOM SHELF DIMENSIONS
9.7"W x 16.5" D

COMPACT 40 WITH 250mL FLASKS




CALCULATE:
Approximately 11
Flask Bottles (250mL) size will
fit on Bottom shelf.




BOTTOM SHELF DIMENSIONS
9.7" W x 16.5" D

COMPACT 40 WITH 500mL FLASKS



CALCULATE:

Approximately 8
Flask Bottles (500mL) size will
fit on Bottom shelf.



BOTTOM SHELF DIMENSIONS

9.7" W x 16.5" D

COMPACT 40 WITH TEST TUBE RACKS

Typical Rack Size
Holds 60 Tubes



Tube size
18 - 20 mm

3.25"
82.6 mm


CALCULATE:

Approximately 2 Test
Tube Racks will fit on
Bottom shelf.

9.625"
244.5mm

6.125"
155.6mm


VESSEL




BOTTOM SHELF DIMENSIONS

9.7" W x 16.5" D

COMPACT 40 WITH AUTOCLAVE WASTE BAG



CALCULATE:
Approximately 1 Small
Waste Bag and Tray will
fit on Bottom shelf.



BOTTOM SHELF DIMENSIONS
9.7" W x 16.5" D