

H60 OUTSIDE DIMENSIONS

TOP SHELF  
DIMENSIONS

13.6" W x 24.4" D  
with 6.9" clearance

BOTTOM SHELF  
DIMENSIONS

9.7" W x 24.4" D  
with 11.6" clearance  
( upper shelf removed )


32.1"  
815mm

20.5"  
520mm


25.2"  
640mm

20.5"  
520mm

DOOR SWING




## H60 WITH 1 LITER FLASKS



CALCULATE:



Approximately 7  
Flask Bottles ( 1 L ) size will  
fit on Bottom shelf.



BOTTOM SHELF DIMENSIONS


9.7" W x 24.4" D

## H60 WITH 2 LITER FLASKS




BOTTOM SHELF DIMENSIONS  
9.7" W x 24.4" D

## H60 WITH 125mL FLASKS




CALCULATE:  
Approximately 63  
Flask Bottles ( 125mL ) size will  
fit on Top and Bottom shelves.




TOP SHELF DIMENSIONS  
13.6"W x 24.4"D  
WITH 6.9" CLEARANCE

BOTTOM SHELF DIMENSIONS  
9.7"W x 24.4" D

## H60 WITH 250mL FLASKS




CALCULATE:  
Approximately 17  
Flask Bottles ( 250mL ) size will  
fit on Bottom shelf.




BOTTOM SHELF DIMENSIONS  
9.7" W x 24.4" D

## H60 WITH 500mL FLASKS



CALCULATE:

Approximately 12  
Flask Bottles (500mL ) size will  
fit on Bottom shelf.



BOTTOM SHELF DIMENSIONS

9.7" W x 24.4" D

## H60 WITH TEST TUBE RACKS

Typical Rack Size  
Holds 60 Tubes



Tube size  
18 - 20 mm


3.25"  
82.6 mm

CALCULATE:

Approximately 4 Test  
Tube Racks will fit on  
Bottom shelf.

9.625"  
244.5mm


6.125"  
155.6mm




BOTTOM SHELF DIMENSIONS

9.7" W x 24.4" D

## H60 WITH AUTOCLAVE WASTE BAG



CALCULATE:  
Approximately 1 Small  
Waste Bag and Tray will  
fit on Bottom shelf.



BOTTOM SHELF DIMENSIONS  
9.7" W x 24.4" D