











Illustrated Spares Price List



2/7/2022

Part Number	Description	List Price
ADT/TEE/015	Pipe Connector - Tee 15mm Compression	\$13.90 
AFL/ACS/EXH	Containment Filter Element - External <i>For Filter Housing</i>	\$695.40 
BFI/ACS/015	Tank Connector 15mm Complete Plated	\$23.40 
BFI/ACS/022	Tank Connector 22mm Complete Plated	\$32.10 
BFI/ACS/028	Tank Connector 28mm Finish -	\$38.50 
CAS/MID/000	Castor Wheel Set - Compact C60 model <i>Set of 4</i>	\$44.50 
CON/230/24A	Contactors - 50amp	\$149.80 
CON/240/09A	Contactors - 9 Amp 240V	\$73.60 
CON/240/23A	Contactors - 23 Amp 240V	\$126.50 
CON/240/45A	Contactors - 100 Amp 240V	\$280.30 

CON/OLD/004 Overload Relay (2.4 - 4.0 Amp range)

\$65.10



CON/OLD/69A Overload Relay (6-9 Amp range)

\$115.70




CON/OLD/OL6 Overload Relay (4-6 Amp range)

\$54.10



FAN/ACS/080 Instrument Cooling Fan - All Models

\$58.20



FAN/ACS/120 Cooling Fan -Compact Models 120mm
All C40 and H60 Models and C60 Models after S/No: 3786

\$64.90



FAN/QCS/UUL Cooling Fan Top and Front Loading Models
UL Approved

\$218.80



FIL/SUP/001 Suppressor 0.1 uF

\$12.60



FIL/VDR/002 Varistor (VDR)

\$5.90




FLO/WAT/SWU Float Switch UL - Vac Header Pump Tank

\$54.60




GAS/ACS/500 Door Seal Circular 100-200L Models

\$176.00



GAS/ACS/Q63 Door Seal Circular 320-400L Models

\$190.40



GAS/QCA/350 Door Gasket QCA Compact Series

\$121.90



GAT/Q63/SAF Q63 Safety Gate

\$120.50




GAT/QCS/SAF Safety Gate QCS

\$88.00



HCO/ACS/150 Thermal Fuse 140°C

\$8.90



HND/QCS/YEL Door Handle - QCS & Compact Series

\$64.50



HTR/227/C40 Heater 320 and 400L Models 2.9KW 277V

\$198.00

Please note: Both leads to have RED markers. INCOLOY-800 277V - 2900W 170°C UL (US models only)



HTR/227/HOR Heater EH Models 3.08KW 277V

\$198.00

Please note: Both leads to have RED markers. INCOLOY-800 277V - 3150W 170°C UL Vertical Heater for US models only.



HTR/227/VER Heater EV Models 3.15KW 277V

\$198.00

Please note: Both leads to have RED markers. INCOLOY-800 277V - 3150W 170°C UL Vertical Heater for US models only.



HTR/ACS/HOR Heater - EH Models 3.5KW 230V

\$198.00



HTR/ACS/VER Heater - EV Models 3.5KW 230V

\$198.00



HTR/MAT/001 Heater Mat 350W 230V - Rectangular

\$111.50

UL Approved



HTR/MAT/002 Heater Mat 350W 230V - Square

\$111.50

UL Approved



HTR/MID/C40 Heater Benchtop & 320-400L Models 2.9kW

\$203.90

230V 13A



HTR/MID/C60 Heater Top Loading Compact Models 2.9kW
230V 13A

\$203.90



INL/ACS/FIL Inline Electrical Filter

\$50.30



KEY/SLK/001 Setting Lock Key Switch

\$133.80



KEY/SLK/003 Setting Lock Key ONLY

\$49.80



KEY/STP/455 Key - Emergency Stop

\$28.50



KEY/TLK/002 Thermal Lock Key Switch - TACTROL2

\$69.50



KEY/TLK/004 Thermal Lock Key ONLY - TACTROL2

\$49.80



LOO/QCS/DOR Door Wiring Loom QCA/QCS

\$330.90



LWP/ACS/00C Low Water Probe Type C all Models
All Models

\$229.20



MCB/UL2/32B Circuit Breaker 32Amp 2Pole Type B UL
UL Approved

\$97.50



MCB/UL2/B06 Circuit Breaker 6Amp 2 Pole Type B UL
UL Approved

\$57.50



MCB/UL2/B10 Circuit Breaker 10Amp 2 Pole Type B UL
UL Approved

\$59.20



MCB/UL2/B13 Circuit Breaker 13Amp 2 Pole Type B UL
UL Approved

\$59.20



MCB/UL3/05C Circuit Breaker Type C-3Pole 5Amp UL
UL Approved

\$118.20



MCB/UL3/B25 Circuit Breaker 25Amp 3 Pole Type B UL
UL Approved

\$70.30



MCB/UL3/C25 Circuit Breaker 25Amp 3 Pole Type C UL
UL Approved

\$70.30



MCB/UL3/C40 Circuit Breaker 40Amp 3 Pole Type C UL
UL Approved

\$65.90



MSA/DCL/QCS Microswitch Cross Roller QCS Door Closed

\$79.00



MSA/LID/QCS Microswitch Small
Top and Front Loaders Door/ Canopy/Bench-Top and Compact Door

\$19.30



NIP/ACS/012 Nipple Brass Hexagon 1/2in MxM

\$3.90



NIP/ACS/034 Nipple Brass Hexagon 3/4in MxM

\$6.70



NIP/ACS/100 Nipple Hexagon Brass 1in MxM

\$7.40



PNC/COM/PWF Power Cord Rated 20Amp/250V
L6-20P to C19 Power Cord, 20A, 250V, 12/3 SJT, 8 Foot

\$69.10















PRB/ACS/003 Thermocouple Probe Duplex K Type














\$146.40

















PRB/ACS/004	Thermocouple Probe Wandering K Type <i>Please note: Wandering probes are sold packs of 2</i>	\$334.30	
PRB/ACS/005	Thermocouple Probe Load Simulator K Type	\$75.10	
PRB/ACS/FSS	Fitting Bulkhead Duplex ProbeT/C 1/8 S/S <i>For thermocouples</i>	\$15.20	
PRB/PLU/SOC	Plug & Socket K Type Wandering Probe	\$15.20	
PRB/PRT/003	Probe Duplex PT100 6 x70mm <i>6mm replaces 3mm probe from Serial No 4800. For first change on models before serial No. 4800 a new fitting (PRB/PRT/FIT) must be supplied</i>	\$129.30	
PRB/PRT/004	PRT Probe Wandering <i>Please note: Wandering probes are sold packs of 2</i>	\$334.30	
PRB/PRT/005	PRT Probe Load Simulator	\$93.00	
PRB/PRT/007	Probe Duplex PT100 6mm x 250mm <i>6mm replaces 3mm probe from Serial No 4800. For first change on models before serial No. 4800 a new fitting (PRB/PRT/FIT) must be supplied</i>	\$129.30	
PRB/SST/PRB	Grade 316 Stainless Steel Cable Tie Pk3 <i>For attaching heater protection probes to heaters. Size: 127 x 4.6 mm. Requires specialist fitting tool.</i>	\$184.90	
PSG/ACS/100	Pressure Gauge 100mm (-1 to 3 bar)	\$230.50	
PSG/GAU/063	Pressure Gauge 63mm (-1 to 3 bar) <i>Bench-Top and Compact Models</i>	\$84.10	
PSS/SWI/001	Pressure Switch Range 1.5 to 15 psi	\$167.30	

PSS/SWI/FIT	Fitting Compression Str.1/8BSP Fim x1/8 <i>Pressure Switch Adaptor - 1/8 MPT to Capillary Pipe (Compression joint).</i>	\$11.50	
PST/ACS/001	Pressure Transducer 0 to 4 Bar <i>Autoclave Pressure</i>	\$258.00	
RLY/MID/HTR	Relay Heater Control - Compact Series	\$57.40	
SBG/ACS/015	Silicon Bung 15mm	\$4.60	
SBG/ACS/022	Silicon Bung 22mm	\$7.90	
SHM/COM/SET	Shim Set Bench-Top & Compact Models <i>Contents: 5 of each shim, Hinge and Clamp 1mm and 0.5mm</i>	\$105.70	
SHM/Q63/SET	Shim Set 320-400L Models <i>Contents: 5 of each shim, Hinge and Clamp 1mm and 0.5mm</i>	\$105.70	
SHM/QCS/SET	Shim Set 100-200L Models <i>Contents: 5 of each shim, Hinge and Clamp 1mm and 0.5mm</i>	\$105.70	
SHM/RSC/001	Shim Hinge Block 320-700L Models 1mm <i>Stainless Steel</i>	\$5.30	
SHM/RSC/005	Shim Hinge Block 320-700L Models 0.5mm <i>Stainless Steel</i>	\$5.20	
SOL/MID/001	12V Door Locking Solenoid Kit. 40-85L <i>Replaces 240V Solenoid on Compact Models before Serial Number 4438. Includes 12V Solenoid, Actuators, Bracket, Microswitch and 12V power supply</i>	\$330.90	
SOL/QCS/001	12V Door Locking Solenoid Kit.100-400L <i>Replaces Door Solenoid on QCS and Q63 Models before Serial Number 4438. Includes 12V Solenoid, Actuator and Bracket, If replacing a 240V solenoid a power supply (TAC/POW/002) is required and must be ordered separately.</i>	\$185.50	

SOL/QCS/002	Door Lock Solenoid 12V 40-400L <i>For Top and Front Loading, Bench-Top and Compact Models - Solenoid Only. Serial Nos from 4438.</i>	\$220.10	
STP/ACS/001	Emergency Stop Button <i>Without Contact Blocks</i>	\$122.60	
STP/ACS/002	Emergency Stop Contact Block	\$62.30	
STR/ACS/015	Strainer Brass 15mm	\$79.50	
STR/ACS/022	Strainer Brass 22mm	\$127.00	
STR/ACS/B22	Mesh Basket 3/4" Drain Filters <i>For threaded fittings only</i>	\$40.40	
STR/ACS/WHL	Conical Strainer Large <i>Fits 28mm Tank Connectors</i>	\$45.90	
STR/ACS/WHM	Conical Strainer Medium <i>Fits 22mm Tank Connectors</i>	\$40.40	
STR/ACS/WHs	Conical Strainer Small <i>Fits 15mm Tank Connectors</i>	\$35.70	
TAC/OPT/001	Tactrol Output Board	\$322.00	
TAC/OPT/FUS	Fuse 3.15A 5x20mm Tactrol O/P Board	\$9.80	
TAC/PCB/001	Tactrol Main Board	\$2,262.00	

TAC/PCB/MIN	Tactrol Mini Processor Main Board	\$836.20	
TAC/POW/001	Power Supply Switch Mode	\$359.90	
TAC/POW/002	Power Supply 240V-12V Solenoid Door	\$74.00	
TAC/PRN/000	Tactrol Thermal Printer <i>Includes Thermal Printer, Serial Interface Adaptor and Mounting Plate</i>	\$1,255.20	
TAC/RIB/DIS	Ribbon Cable Display to Main Board Long <i>Models with Tactrol Board Not Mounted in Control Box</i>	\$76.10	
TAC/RIB/DISS	Ribbon Cable Display to Main Board Short <i>Models with Tactrol Board Mounted in Control Box</i>	\$38.30	
TAC/RIB/PRIL	Ribbon Cable Impact Printer Long <i>Models without a control box</i>	\$72.60	
TAC/RIB/PRIS	Ribbon Cable Impact Printer Short <i>Models with a control box</i>	\$40.30	
TAC/RIB/PRN	Ribbon Cable Printer to Main Board <i>Power Door and Double Door Models. Serial Number required - length varies with model</i>	\$56.40	
TAC/RIB/PRNL	Ribbon Cable Thermal Printer Long <i>Models without a control box - Thermal Printer</i>	\$72.60	
TAC/RIB/PRNS	Ribbon Cable Thermal Printer Short <i>Models with a control box</i>	\$40.30	
TAC/USB/001	Tactrol USB Interface Card	\$358.50	
TAC/USB/COV	Tactrol USB Black Dust Cap	\$10.50	

TC2/CON/001	Tactrol 2 Main Control Panel	\$717.20	
TC2/PRG/005	Tactrol 2 Programmer Panel	\$871.00	
THW/ACS/002	Load Simulator Thermowell Compact /QCS	\$132.60	
TST/STG/210	Thermostat 40 to 210C <i>Drain Condensers and Steam Generators</i>	\$55.50	
VAC/LRP/LRG	Vacuum Pump Liquid Ring 3 phase 450-700L <i>Q63 or RSC Models</i>	\$3,026.60	
VAC/LRP/MD2	Vacuum Pump 1Phase 60 Hz UL 100-400L <i>UL Approved US Market only</i>	\$2,537.50	
VAC/LRP/MED	Vacuum Pump Liquid Ring 3 phase 230-400L	\$2,692.10	
VAC/LRP/SM3	Vacuum Pump Liquid Ring 3 phase 100-200L <i>100-200L Models</i>	\$2,225.90	
VCE/ACS/001	Condensate Trap BPT13S 1/2in 15mm	\$272.30	
VFL/ACS/TOR	Valve Torbeck Type Auto Water Fill Tank	\$46.40	
VNR/ACS/015	Valve Non-Return 15mm <i>(1/2in)</i>	\$36.70	
VNR/ACS/022	Valve Non-Return 22mm <i>(3/4in)</i>	\$44.40	

VNR/ACS/028 Valve Non-Return 28mm (1in)	\$50.60	
VPC/ACS/001 Air Purge Valve	\$349.10	
VSC/ASM/015 Safety Valve 15mm, 2.4Bar US Models (UL) <i>Approval acc. to ASME Boiler and Pressure Vessel Code, Section VIII.Div 1 (ASME)</i>	\$384.30	
VSC/SSM/015 Safety Valve 1/2" x1" Stainless ASME <i>Approval acc. to ASME Boiler and Pressure Vessel Code, Section VIII.Div 1 (ASME)</i>	\$395.90	
VSI/ACS/015 Steam Inlet Valve 15mm	\$253.80	
VVC/RSC/028 Valve Solenoid 28mm <i>For serial numbers lower than 4628 please order the conversion kit VVC/RSC/ADP to ensure the correct valve length for direct replacement.</i>	\$329.70	
VVC/RSC/ADP 28mm Solenoid Valve Converter Kit <i>Adapter Kit to fit 28mm Solenoid Valve to Autoclaves before serial number 4628.</i>	\$45.70	
VVC/UUL/015 Valve Solenoid 15mm UL Approved	\$185.80	
VVC/UUL/022 Valve Solenoid 22mm UL Approved	\$199.00	
VVC/UUL/028 Valve Solenoid 28mm UL Approved <i>For serial numbers lower than 4628 please order the conversion kit VVC/RSC/ADP to ensure the correct valve length for direct replacement.</i>	\$337.80	
VVF/ACS/015 Fitting Vent Valve 15mm x 1/2 MIC	\$6.50	

VVF/ACS/022 Fitting Vent Valve 22mm x 3/4 MIC

\$7.30



VVF/ACS/028 Fitting Vent Valve 28mm x 1in MIC

\$11.90



WAS/PAC/015 Fibre Washers 15mm Pack of 10

\$1.60



WAS/PAC/022 Fibre Washers 22mm Pack of 10

\$1.60



WAS/PAC/028 Fibre Washers 28mm Pack of 10

\$1.60



WAS/SET/001 Fibre Washer Set

Complete set of washers, with PTFE tape & jointing compound

\$86.30



ZZZ/RIB/002 Tactrol Printer Ribbon Type 2

\$34.60



ZZZ/ROL/005 Printer Roll (Pk. 5)

\$34.60



ZZZ/ROL/007 Thermal Printer Roll (Pk 5)

\$41.60



ZZZ/SYP/LGE Siphon Pump Large

Long Siphon Pump for Top Loading Autoclaves

\$74.00



ZZZ/SYP/SML Siphon Pump Small

Short Siphon Pump for Front Loading Autoclaves

\$44.00



ZZZ/TAP/001 Tape PTFE 12Mx12mm roll

\$10.30



ZZZ/TCG/001 Grease Copper Tube

\$41.50



ZZZ/TSG/001 Grease Silicone Tube 100g

\$48.00

