

Autoclave Model No.	Supply	Autoclave Amps	ACME Model	VA	Branch MCB
PNA/QCS/C40	phase to phase or single phase	14	T-1-81051	500	15A
PNA/QCS/C60					
PNA/QCS/H60					
PNA/QVA/C40					
PNA/QVA/C60					
PNA/QVA/H60					
PNA/QVA/C40M					
PNA/QVA/C60M					
PNA/QVA/H60M					
PNA/QCS/EV100					
PNA/QCS/EV150					
PNA/QCS/EH100					
PNA/QCS/EH150					
PNA/QCS/EH150D					
PNA/QCS/EH200					
PNA/QCS/EH200D					
PNA/QVA/EV100					
PNA/QVA/EV100M					
PNA/QVA/EV150					
PNA/QVA/EV150M					
PNA/QVA/EH100					

Autoclave Model Number	Supply	Autoclave Amps	ACME Model	VA	Branch MCB
PNA/QVA/EH100M	3 phase	28	T-1-81052 x 2	750 x 2	30A
PNA/QVA/EH150					
PNA/QVA/EH150M					
PNA/QVA/EH150D					
PNA/QVA/EH150DM					
PNA/QVA/EH200					
PNA/QVA/EH200M					
PNA/QVA/EH200D					
PNA/QVA/EH200DM					
PNA/QCS/SV100	3 phase	3	T-1-81051 x 2	500 x 2	5A
PNA/QCS/SV150					10A
PNA/QCS/SH100					5A
PNA/QCS/SH150					10A
PNA/QCS/SH150D					5A
PNA/QCS/SH200					10A
PNA/QCS/SH200D					5A
PNA/QCS/SV100					phase to phase or single phase
PNA/QCS/SV150					
PNA/QCS/SH100	10.5	15A			
PNA/QCS/SH150					
PNA/QCS/SH150D					
PNA/QCS/SH200					
PNA/QCS/SH200D					
PNA/QVA/SV100	3 phase	11	T-1-81051 x 2	500 x 2	15A
PNA/QVA/SV100M		15			20A
PNA/QVA/SV150		11			15A
PNA/QVA/SV150M		15			20A

Autoclave Model Number	Supply	Autoclave Amps	ACME Model	VA	Branch MCB
PNA/QVA/SH100	3 Phase	11	T-1-81051 x 2	500 x 2	15A
PNA/QVA/SH100M		15			20A
PNA/QVA/SH150		11			15A
PNA/QVA/SH150M		15			20A
PNA/QVA/SH150D		11			15A
PNA/QVA/SH150DM		15			20A
PNA/QVA/SH200		11			15A
PNA/QVA/SH200M		17.5			20A
PNA/QVA/SH200D		11			15A
PNA/QVA/SH200DM		17.5			20A
PNA/QCS/EH320		3 Phase			45
PNA/QCS/EH400					
PNA/QCS/EH320	phase to phase or single phase	73	T-1-11685	2000	80A
PNA/QCS/EH400					
PNA/QVA/EH320	3 Phase	45	T-1-11684 x 2	1500 x 2	50A
PNA/QVA/EH320M					
PNA/QVA/EH400					
PNA/QVA/EH400M					
PNA/QVA/EH320	phase to phase or single phase	75	T-1-11685	2000	80A
PNA/QVA/EH320M					
PNA/QVA/EH400					
PNA/QVA/EH400M					
PNA/QCS/SH320	phase to phase or single phase	5	T-1-81051	500	8A
PNA/QCS/SH500					
PNA/QVA/EH320	phase to phase or single phase	16	T-1-81051	500	20A
PNA/QVA/EH400					
PNA/QVA/EH320M		27	T-1-81052	750	30A
PNA/QVA/EH400M					

Note: Circuit breakers protecting the transformer to be set at not more than 125% of the rated full-load input current of the transformer in accordance with The National Electrical Code 450-4.

